



How to install the BarTender UltraLite software

1. Locate the software installation file, BT2016_UL, on this flash drive and start the installation.
2. Follow the instructions on screen to install the software.

How to install the DiTi™ 203iD driver

1. To install the required printer driver, locate the “2021.2_M-4” file in the “Drivers” folder to start the installation process.
2. Follow the instructions on screen until the installation process is completed successfully.

Please contact us for technical support at: support@ga-international.com

or call us toll free at **1-800-518-0364 (US and Canada)**.

+1 (450) 473-9420 (International).

ga-international.com

Labtag.com