



Autoclave Labels & Tapes

Labtag.com/autoclave

Steri-LazrTAG™

Autoclave-Resistant Labels for Laser Printers

Temperature range -40°C to +150°C (-40°F to +302°F)

These autoclave-resistant laser labels withstand multiple autoclave cycles and are ideal for identifying materials that will undergo repeated sterilization protocols. They are available in permanent or removable, transparent, and self-laminating, wrap-around formats for extra durability.

Permanent				
SKU#	Inches (W x H)	mm (W x H)	Labels/Pack	Color Available
AKA-12	0.94" x 0.5"	23.9 x 12.7	2,688	Ⓜ
AKA-13	1.75" x 0.5"	44.5 x 12.7	1,280	Ⓜ
AKA-63	2" x 1.125"	50.8 x 28.6	432	Ⓜ
AKA-3	2.63" x 1"	66.8 x 25.4	480	Ⓜ
AKA-30	3" x 2.375"	76.2 x 60.3	128	Ⓜ
AKA-28	8.5" x 11"	215.9 x 279.4	16	Ⓜ
Removable				
SKU#	Inches (W x H)	mm (W x H)	Labels/Pack	Colors Available
ALTR-13	1.75" x 0.5"	44.5 x 12.7	1,280	Ⓜ YE GA BL OR RE LA
ALTR-63	2" x 1.125"	50.8 x 28.6	432	Ⓜ
ALTR-3	2.63" x 1"	66.8 x 25.4	480	Ⓜ
ALTR-116	3" x 0.625"	76.2 x 15.9	512	Ⓜ
ALTR-39	3.5" x 1"	88.9 x 25.4	320	Ⓜ OR



Steri-JetTAG™

Autoclave-Resistant Labels for Inkjet Printers

Temperature range -40°C to +150°C (-40°F to +302°F)

These high-temperature resistant permanent inkjet labels are also resistant to spraying with ethanol, isopropanol, and other alcohols and allow full-color text and images to be printed.

SKU#	Inches (W x H)	mm (W x H)	Labels/Pack
INKA-13	1.75" x 0.5"	44.5 x 12.7	1,280
INKA-63	2" x 1.125"	50.8 x 28.6	432
INKA-106	3" x 2"	76.2 x 50.8	160
INKA-61	3.5" x 2"	88.9 x 50.8	160



Print-N-Shield™

Self-Laminating Labels for Autoclave

Temperature range -60°C to +150°C (-76°F to +302°F)

Self-laminating, Print-N-Shield labels come with a clear laminating tab that protects the label and its printout from water, detergents, abrasion, harsh chemicals, and storage in a wide range of temperatures. Perfect for labeling equipment and material/containers that will undergo routine cleaning or steam autoclaving.

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll	Color Available
SLPA-534NPC1-0.5WH	0.75" x 2.25"	19.1 x 57.2	1"	500	Ⓜ
SLPA-534NPC3-1WH	0.75" x 2.25"	19.1 x 57.2	3"	1000	Ⓜ
SLPA-284C3-1WH	1" x 2.125"	25.4 x 54	3"	1000	Ⓜ



Steri-ThermoTAG™

Autoclave-Resistant Thermal-Transfer Labels

Temperature range -60°C to +150°C (-76°F to +302°F)

Thermal-transfer autoclavable labels endure various sterilization processes, from steam autoclaves and dry heat ovens to ethylene oxide gas and gamma irradiation. Blackout Steri-ThermoTAG are also available for concealing sensitive information on containers that need to be autoclaved prior to disposal.



Permanent

SKU# (White)	SKU# (Colors)	Inches (W x H)	mm (W x H)	Core	Labels/Roll	Colors Available
AUTT-75NPC1-2	AUTC-75NPC1-2	1" x 0.6875"	25.4 x 17.5	1"	2,000	WH GA OR
AUTT-75NPC3-2	AUTC-75NPC3-2	1" x 0.6875"	25.4 x 17.5	3"	2,000	WH GA OR
AUTT-19C1-1WH	-	1.3" x 0.625"	33 x 15.9	1"	1,000	WH
AUTT-19C3-1WH	-	1.3" x 0.625"	33 x 15.9	3"	1,000	WH
AUTT-84C1-0.25	AUTC-84C1-0.25	2" x 4"	50.8 x 101.6	1"	250	WH YE GA PI BL OR RE
AUTT-84C3-1	AUTC-84C3-1	2" x 4"	50.8 x 101.6	3"	1,000	WH YE GA PI BL OR RE
AUTT-28C1-2WH	-	2" x 1"	50.8 x 25.4	1"	2,000	WH
AUTT-28C3-2WH	-	2" x 1"	50.8 x 25.4	3"	2,000	WH
AUTT-200C1-1WH	-	2.5" x 1.25"	63.5 x 31.8	1"	1,000	WH
AUTT-200C3-2WH	-	2.5" x 1.25"	63.5 x 31.8	3"	2,000	WH
AUTT-62C1-0.5WH	-	4" x 2"	101.6 x 50.8	1"	500	WH
AUTT-62C3-1WH	-	4" x 2"	101.6 x 50.8	3"	1,000	WH

Removable

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll	Color Available
AUTR-61NPNOTC1-4WH	0.354" circle	9	1"	4,000	WH
AUTR-61NPNOTC3-8WH	0.354" circle	9	3"	8,000	WH
AUTR-9C1-1WH	1.25" x 0.875"	31.8 x 22.2	1"	1,000	WH
AUTR-9C3-3WH	1.25" x 0.875"	31.8 x 22.2	3"	3,000	WH
AUTR-28C1-1WH	2" x 1"	50.8 x 25.4	1"	1,000	WH
AUTR-28C3-2WH	2" x 1"	50.8 x 25.4	3"	2,000	WH
AUTR-11C1-2WH	2.64" x 0.276"	67.1 x 7	1"	2,000	WH
AUTR-11C3-4WH	2.64" x 0.276"	67.1 x 7	3"	4,000	WH

Blackout

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll	Color Available
BOA-19C1-1SI	1.3" x 0.625"	33 x 15.9	1"	1,000	SI
BOA-19C3-3SI	1.3" x 0.625"	33 x 15.9	3"	3,000	SI
BOA-28C1-1SI	2" x 1"	50.8 x 25.4	1"	1,000	SI
BOA-28C3-2SI	2" x 1"	50.8 x 25.4	3"	2,000	SI

Autoclave Piggyback Labels

Autoclave-Resistant Thermal-Transfer Piggyback Labels

Temperature range -60°C to +150°C (-76°F to +302°F)

These thermal-transfer printable labels have a label-on-label multi-layer design that provides a versatile identification solution for various containers and surfaces. The smaller nested labels can be used to identify different surfaces or containers as needed, providing additional tracking options.



SKU#	Total set size (W x H) Inches / mm	Main Label (1) (W x H) Inches / mm	Secondary Labels (4) (W x H) Inches / mm	Core	Labels/Roll
PAUTR-2C1-0.5WH	4" x 2" / 101.6 x 50.8	2.5" x 1.75" / 63.5 x 44.5	1.25" x 0.25" / 31.8 x 6.4	1"	500
PAUTR-2C3-0.5WH	4" x 2" / 101.6 x 50.8	2.5" x 1.75" / 63.5 x 44.5	1.25" x 0.25" / 31.8 x 6.4	3"	500
PAUTR-6WRC1-0.5WH	4" x 2" / 101.6 x 50.8	2.5" x 1.75" / 63.5 x 44.5	1.25" x 0.375" / 31.8 x 9.5	1"	500
PAUTR-6WRC3-0.5WH	4" x 2" / 101.6 x 50.8	2.5" x 1.75" / 63.5 x 44.5	1.25" x 0.375" / 31.8 x 9.5	3"	500

Steam Sterilization Indicator Labels

Autoclave-Resistant Labels with Indicator Strip

Temperature range -60°C to +150°C (-76°F to +302°F)

These permanent autoclave-resistant labels are coated with an indicator strip that will noticeably change color (white to charcoal gray) following a steam autoclave sterilization cycle. This allows material to be clearly and accurately identified, while also providing a visual indicator of autoclaving.



SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll
SAUP-56C1-0.5WC	3" x 2"	76.2 x 50.8	1"	500
SAUP-56C3-1WC	3" x 2"	76.2 x 50.8	3"	1,000

Wrap-Around Cryo-LazrTAG™

Cryogenic & Autoclave Resistant Self-Laminating Laser Labels

Temperature range -196°C to +150°C (-321°F to +302°F)

These self-laminating cryogenic laser labels are perfect for use in extreme temperatures and in various harsh environments. The clear laminating part of the labels protects the printed information from autoclaving, boiling water, cryo storage, abrasion, and chemical exposure.

SKU#	Inches (W x H)	mm (W x H)	Labels/Pack	Colors Available
DFSL-13WH	1" x 0.5" + 0.75"	25.4 x 12.7 + 19.1	3,200	ⓂⓂ
DFSL-62WH	1" x 0.6" + 0.625"	25.4 x 15.2 + 15.9	2,560	ⓂⓂ
DFSL-65	1" x 0.6" + 1.6"	25.4 x 15.2 + 40.6	1,920	ⓂⓂ ⓂⓂ ⓂⓂ ⓂⓂ
DFSL-75WH	1" x 0.625" + 2"	25.4 x 15.9 + 50.8	1,280	ⓂⓂ
DFSL-77WH	2" x 0.625" + 2"	50.8 x 15.9 + 50.8	1,280	ⓂⓂ



Transparent Permanent Cryo-LazrTAG™

Cryogenic & Autoclave Resistant Laser Labels

Temperature range -196°C to +150°C (-321°F to +302°F)

Our transparent Cryo-LazrTAG cryogenic laser labels will also withstand sterilization cycles in steam autoclaves, as well as gamma irradiation, up to 50 KGy. Their clear design ensures the contents of labeled containers remain clearly visible, and can be inscribed with permanent markers.

SKU#	SKU# (with Color Tab)	Inches (W x H)	mm (W x H)	Labels/Pack	Color Tab
DFLT-121	-	0.394" x 0.315"	10 x 8	7,920	Transparent
DFLT-62	-	1.625" x 0.6"	41.3 x 15.2	1,024	Transparent
DFLT-13	DFLT-13B	1.75" x 0.5"	44.5 x 12.7	1,280	ⓂⓂⓂ ⓂⓂⓂ ⓂⓂⓂ ⓂⓂⓂ
DFLT-65	-	2.6" x 0.6"	66 x 15.2	768	Transparent
DFLT-75	-	3" x 0.625"	76.2 x 15.9	512	Transparent
DFLT-77	-	4" x 0.625"	101.6 x 15.9	512	Transparent



Cryo-OmniTAG™

Transparent Cryogenic & Autoclave-Resistant Labels

Temperature range -196°C to +150°C (-321°F to +302°F)

These versatile thermal-transfer labels are well suited for a wide range of applications, from cryogenic storage in liquid nitrogen to autoclaving and dry-heat sterilization. Their transparent design ensures they will not conceal the contents of labeled containers, perfect for cryo vials, freezer boxes, glass bottles, and more.

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll
GANA-52NOTC1-2	0.94" x 0.5" + 0.437" circle	23.9 x 12.7 + 11.1	1"	set 2,000
GANA-52NOTC3-2	0.94" x 0.5" + 0.437" circle	23.9 x 12.7 + 11.1	3"	set 2,000
GANA-158NOTC1-2	1.25" x 0.5" + 0.437" circle	31.8 x 12.7 + 11.1	1"	set 2,000
GANA-158NOTC3-2	1.25" x 0.5" + 0.437" circle	31.8 x 12.7 + 11.1	3"	set 2,000
GANA-9NOTC1-1	1.25" x 0.875"	31.8 x 22.2	1"	1,000
GANA-9NOTC3-2	1.25" x 0.875"	31.8 x 22.2	3"	2,000
GANA-19NOTC1-1	1.3" x 0.625"	33 x 15.9	1"	1,000
GANA-19NOTC3-2	1.3" x 0.625"	33 x 15.9	3"	2,000
GANA-10NOTC1-1	1.5" x 0.75"	38.1 x 19.1	1"	1,000
GANA-10NOTC3-2	1.5" x 0.75"	38.1 x 19.1	3"	2,000
GANA-107NOTC1-1	1.5" x 0.875"	38.1 x 22.2	1"	1,000
GANA-107NOTC3-3	1.5" x 0.875"	38.1 x 22.2	3"	3,000
GANA-28NOTC1-1	2" x 1"	50.8 x 25.4	1"	1,000
GANA-28NOTC3-2	2" x 1"	50.8 x 25.4	3"	2,000
GANA-210NOTC1-1	2.72" x 1"	69 x 25.4	1"	1,000
GANA-210NOTC3-1	2.72" x 1"	69 x 25.4	3"	1,000



Cryo-WrapTAG™

Wrap-Around Cryogenic & Autoclave-Resistant Labels

Temperature range -196°C to +150°C (-321°F to +302°F)

These self-laminating, wrap-around labels have a clear laminate that is applied over the printed portion of the label, protecting it against harsh lab environments and extreme temperatures. This includes storage in liquid nitrogen, steam autoclaving, as well as abrasion, fading, and exposure to harsh chemicals.

SKU#	Inches (W x H + wrap)	mm (W x H + wrap)	Core	Labels/Roll	Colors Available
CATT-318SBNOTC1-1 WH	1" x 0.25" + 0.75"	25.4 x 6.4 + 19.1	1"	1,000	WH
CATT-318SBNOTC3-2WH	1" x 0.25" + 0.75"	25.4 x 6.4 + 19.1	3"	2,000	WH
CATT-312NOTC1-1WH	1" x 0.625" + 1.375"	25.4 x 15.9 + 35	1"	1,000	WH
CATT-312NOTC3-2 WH	1" x 0.625" + 1.375"	25.4 x 15.9 + 35	3"	2,000	WH
CATT-311NOTC1-2WH	1" x 0.625" + 1"	25.4 x 15.9 + 25.4	1"	2,000	WH
CATT-311NOTC3-4WH	1" x 0.625" + 1"	25.4 x 15.9 + 25.4	3"	4,000	WH
CATT-303NOTC1-1WH	1.125" x 1" + 1.595"	28.6 x 25.4 + 40.5	1"	1,000	WH
CATT-303NOTC3-1WH	1.125" x 1" + 1.595"	28.6 x 25.4 + 40.5	3"	1,000	WH
CATT-308NOTC1-1WH	1.34" x 1" + 1.38"	34 x 25.4 + 35	1"	1,000	WH
CATT-308NOTC3-1WH	1.34" x 1" + 1.38"	34 x 25.4 + 35	3"	1,000	WH
CATT-309NOTC1-1	1" x 1.18" + 1.69"	25.4 x 30 + 42.9	1"	1,000	WH YE GA PI
CATT-309NOTC3-1	1" x 1.18" + 1.69"	25.4 x 30 + 42.9	3"	1,000	WH YE GA PI
CATT-312NOTC1-1WH	1" x 0.625" + 1.375"	25.4 x 15.9 + 35	1"	1,000	WH
CATT-312NOTC3-2WH	1" x 0.625" + 1.375"	25.4 x 15.9 + 35	3"	2,000	WH



Steam Sterilization Indicator Tape

Steam sterilization indicator tape is designed to identify containers that will undergo sterilization in steam autoclaves. This autoclave tape comes with nearly invisible marking strips that turn black when put through a sterilization cycle, which are used to confirm that the labeled material has been autoclaved. Accepts writing with permanent markers.

SKU#	Inches x Feet	mm x m	Core
STRAT-18WH	0.7" x 180'	18mm x 54m	3"
STRAT-25WH	1" x 180'	25.4mm x 54m	3"



#STRAT-25

Before:



After:



Autoclave-Resistant Tape

Temperature range -60°C to +150°C (-76°F to +302°F)

Thermal-transfer printable sterilization tape is designed for use with steam autoclaves and dry heat sterilization. This autoclave-resistant tape is ideally suited for identifying surgical instruments and trays, double-blister packs, as well as plasticware and glassware that will undergo sterilization. Supplied as a continuous roll to provide you with the precise length required.

SKU#	Inches x Feet	mm x m	Core
TAUT-25C1-50WH	1" x 50'	25.4mm x 15m	1"
TAUT-25C3-100WH	1" x 100'	25.4mm x 30m	3"



#TAUT-25

Removable SteriTape™ – Autoclave-Resistant Tape

Temperature Range -80°C to +150°C (-112°F to +302°F)

Removable, autoclave-resistant tape is ideal for identifying reusable sterilization containers intended for use in steam autoclaves, such as Genesis® containers, as well as single-use sterilization bags and pouches. Its removable adhesive ensures it will not leave any adhesive residue upon removal.

SKU#	Inches x Feet	mm x m	Core
TARG-25C3-100WH	1" x 100'	25.4mm x 30m	3"



#TARG-25



+1-450-973-9420
 cx@ga-international.com
 Labtag.com

USA
 P.O. Box 2633,
 Champlain, NY
 12919, United States
 1-800-518-0364

CANADA
 4455 Rue Louis-B.-Mayer,
 Laval QC,
 H7P 6B5, Canada
 1-800-518-0364

NETHERLANDS
 Runweg 28A,
 5258 BN Berlicum
 Netherlands
 +31 (0) 73-7370185

UK
 10 John Street,
 London WC1N 2EB
 United Kingdom
 +44 (20) 3769 5683