

Autoclave Piggyback Labels

For Thermal-Transfer Printers

Temperature range from -60°C to +150°C

- Removable autoclave resistant piggyback labels ideal for metal, plastic, and glass surfaces including surgical instruments, glassware, trays, blister packs, and other.
- Items intended for use in steam autoclaves, ovens, and gamma irradiators. Our durable autoclave thermal-transfer labels will not leave any adhesive residue after removal.
- They are weather & water-proof, alcohol-resistant, and have been tested to endure up to 20 autoclave cycles (up to +132°C / +269.6°F).
- Compatible with all major brands of thermal transfer printers.



SHOP NOW

SKU#	Total Size Inch (W x H)/mm (W x H)	1 Primary Label Inch (W x H)/mm (W x H)	4 Sub-Labels Inch (W x H)/mm (W x H)	Core	Labels/ Roll
PAUTR-2C1-0.5WH	4" x 2" 101.6 x 50.8mm	2.5" x 1.75" 63.5 x 44.5mm	1.25" x 0.25" 31.8 x 6.4mm	1"	500
PAUTR-2C3-0.5WH	4" x 2" 101.6 x 50.8mm	2.5" x 1.75" 63.5 x 44.5mm	1.25" x 0.25" 31.8 x 6.4mm	3"	500
PAUTR-6C1-0.5WH	4" x 2" 101.6 x 50.8mm	2.5" x 1.75" 63.5 x 44.5mm	1.25" x 0.375" 31.8 x 9.5mm	1"	500
PAUTR-6C3-0.5WH	4" x 2" 101.6 x 50.8mm	2.5" x 1.75" 63.5 x 44.5mm	1.25" x 0.375" 31.8 x 9.5mm	3"	500



After printing, peel the label assembly from the main support liner



Affix in your lab notebook or wherever you keep your records



Peel the primary label from the secondary liner



Affix the primary label on the main container or surface being identified



Peel any secondary labels that may be required



Apply the secondary labels to linked containers or documents for tracking purposes