

CalTAG™

Equipment Maintenance & Calibration Labels

Our self-laminating waterproof calibration labels are designed to reliably record the inspection and calibration data on equipment, instruments and containers used in various industries. These labels can be written on with ball-point pens, gel pens, or permanent markers. The transparent tab will securely laminate the data and protect it against tampering, abrasion, and chemical exposure. Available with permanent and removable adhesives.



Step 1: Write



Step 2: Peel



Step 3: Seal



Step 4: Peel



Step 5: Apply

SKU#	Inches (W x H)	mm (W x H)	Printed Text (Title)	Labels/Roll	Adhesive
CALA-001-1R	0.75" x 1.5"	19.1 x 38.1	"CALIBRATION"	1,000	Removable
CALA-001-1P	0.75" x 1.5"	19.1 x 38.1	"CALIBRATION"	1,000	Permanent
CALA-002-1P	1" x 2.125"	25.4 x 54	"CALIBRATION"	1,000	Permanent
CALA-002-1R	1" x 2.125"	25.4 x 54	"CALIBRATION"	1,000	Removable
CALA-005-0.25R	3.5" x 2"	88.9 x 50.8	"SOLUTION"	250	Removable
CALA-006-0.25R	3.5" x 2"	88.9 x 50.8	"DO NOT USE"	250	Removable
CAL-007-1R	1" x 2"	25.4 x 50.8	"OPENED"	1,000	Removable
CALA-008-1R	1" x 2.125"	25.4 x 54	"PREVENTIVE MAINTENANCE"	1,000	Removable
CALA-008-1R-TBL	1" x 2.125"	25.4 x 54	"PREVENTIVE MAINTENANCE"	1,000	Removable
CAL-009-0.25R	3.5" x 2"	88.9 x 50.8	Without Title	250	Removable
CAL-010-1R	1" x 2.125"	25.4 x 54	Without Title	1,000	Removable
CAL-011-0.5R	1.5" x 2.5"	38.1 x 63.5	"OPENED"	500	Removable
CAL-012-0.25R	2.5" x 4"	63.5 x 101.6	"Solution/Reagent"	250	Removable
CALA-014-1R	3.25" x 1"	82.6 x 25.4	"MONTHLY PREVENTIVE MAINTENANCE"	1,000	Removable



Shop Online at [Labtag.com/calibration-labels](https://www.labtag.com/calibration-labels)



#CALA-001

CALIBRATION	
I.D.	
DATE	BY
DUE	

#CALA-008

PREVENTIVE MAINTENANCE	
I.D.	
DATE	BY
DUE	

#CAL-010

Equipment Type:	_____
Serial #	_____
Company Name:	_____
Calibration Date:	_____ dd-mmm-yyyy
Date of expiry:	_____ dd-mmm-yyyy
Initials:	_____ Date dd-mmm-yyyy: _____

#CALA-002

CALIBRATION	
I.D.	
DATE	BY
DUE	

#CALA-008-1R-TBL

PREVENTIVE MAINTENANCE	
I.D.	
DATE	BY
DUE	

#CALA-006

DO NOT USE	
<input type="checkbox"/> Pending Qualification	
<input type="checkbox"/> Pending Repair	
<input type="checkbox"/> Other: _____	
Initial: _____ Date: _____	

#CAL-007

OPENED	
Date: _____ By: _____	
Assigned Expiration:	

#CAL-011

OPENED	
By: _____ Date: _____	
Assigned Expiration:	

#CALA-005

SOLUTION	
Project #:	
Name:	
Lot #:	
Conc.:	pH:
Volume:	
Prepared By:	Date:
Storage Temp:	
EXPIRATION DATE:	

#CAL-009

Project #:	
Name:	
Lot #:	
Conc.:	pH:
Volume:	
Date Received:	By:
Storage Temp:	
EXPIRATION DATE:	

#CAL-012

Solution/Reagent	
Description: _____	
Study/Sponsor #: _____	
Batch/Lot#: _____	
Composition/Concentration: _____	
Prepared By/Date: _____	
Storage: _____	EXPIRATION DATE:

#CALA-014

MONTHLY PREVENTIVE MAINTENANCE	
ID #: _____	
DONE BY: _____	DATE: _____
NEXT MAINT DATE: _____	