

ChemicoTAG™

Solvent-Resistant Color Labels for Containers

Temperature range -20°C to +98°C (-4°F to +208°F)

- Withstand direct exposure to harsh chemicals and solvents
- Durable adhesive ensures labels will not wrinkle or peel
- Available in a variety of colors to increase and simplify identification methods
- Compatible with most thermal-transfer printers



Permanent

These permanent thermal-transfer color labels are perfect for labeling glass and plastic containers that will be used to store harsh chemicals and solvents or will come into contact with them.

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll	Colors
AUC-132C1-0.5	2" circle	50.8 mm circle	C1	500	WH YE BL RE
AUC-85C1-1WH	2" x 0.625"	50.8 x 15.9 mm	C1	1,000	WH
AUC-28C1-1	2" x 1"	50.8 x 25.4 mm	C1	1,000	WH GA
AUC-84C1-0.3WH	2" x 4"	50.8 x 101.6 mm	C1	300	WH
AUC-245C1-1WH	2" x 1.25"	50.8 x 31.8 mm	C1	1,000	WH
AUC-86C1-1	3" x 1.25"	76.2 x 31.8 mm	C1	1,000	WH YE GA RE
AUC-56C1-0.5	3" x 2"	76.2 x 50.8 mm	C1	500	WH YE

Also available on 3" core for industrial printers.



Removable

These removable thermal-transfer labels have a glove-friendly adhesive that ensures the labels can be removed without leaving any adhesive residue behind, ideal for labeling glass, polypropylene, and metal surfaces.

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll	Colors
AUBR-132C1-0.5	2" circle	50.8 mm circle	C1	500	WH YE GA BL RE
AUBR-28C1-1	2" x 1"	50.8 x 25.4 mm	C1	1,000	WH YE GA BL RE PI LA
AUBR-86C1-1	3" x 1.25"	76.2 x 31.8 mm	C1	1,000	WH YE GA RE
AUBR-186C1-0.25	3" x 5"	76.2 x 127 mm	C1	250	WH BL RE
AUBR-62C1-1	4" x 2"	101.6 x 50.8 mm	C1	1,000	WH YE RE LA
AUBR-122SBC1-0.3	4" x 4"	101.6 x 101.6 mm	C1	300	WH YE

Also available on 3" core for industrial printers.

