

Cryo DTermo™

Cryogenic Labels for DYMO® LabelWriter™ 450-Series

Temperature range from -196°C to +70°C




- Cryo DTermo™ labels for long term cryo-preservation in liquid phase and vapor phase liquid nitrogen (-196°C/-321°F)
- Waterproof, can be immersed in liquid nitrogen directly following application, and withstand multiple freeze-thaw cycles
- Able to conform to microplates, freezer boxes, plastic bags, microtubes, microcentrifuge tubes, and cryo vials (0.5 - 5 mL)
- Labels come with downloadable templates to format and print with the DYMO Label Software
- Offered in multiple colors for easy color-coding and customizable with logos, graphics, barcodes and more
- Available in multiple sizes and shape options (rectangle, circle, or rectangle-and-circle set)


Shop on
[Labtag.com/cryodtermo](https://labtag.com/cryodtermo)



Cryo DTermo™ DYMO-Compatible Cryogenic Labels



Item #	Label Size (W x H)	Labels Across	Qty/Roll	Color
EF1F-030	0.35" (9mm) circle	2	2,000	




Item #	Label Size (W x H)	Labels Across	Qty/Roll	Color
EF1F-040	0.5" x 1" (12.7mm x 25.4mm)	2	1,000	



Item #	Label Size (W x H)	Labels Across	Qty/Roll	Color
EF1F-073	0.7" x 1.49" (17.8mm x 37.9mm)	1	500	



Item #	Label Size (W x H)	Labels Across	Qty/Roll	Color
EF1F-070	1" x 1" (25.4mm x 25.4mm)	1	750	



Item #	Label Size (W x H)	Labels Across	Qty/Roll	Color
EG1F-081	1.04" x 0.5" + 0.375" circle (26.4mm x 12.7mm + 9.5mm circle)	1 set	500	

