

Cryo-JetTAG™

Cryogenic Inkjet Labels

- Ideal for prolonged storage in low temperature conditions, including in liquid nitrogen (-196°C) and lab freezers (-80°C)
- Printout has excellent resistance to harsh chemicals such as alcohol, acetone, and xylene
- Available in both sheet and roll formats, for printing with desktop and roll inkjet printers
- Perfect for printing color logos, graphics, text, barcodes, serial numbers, and other types of information



Roll Format

Temperature range -196°C to +121°C (-321°F to +250°F)

Various models of roll format inkjet printers can be used to print these inkjet labels; They can be printed with eye-catching color logos, graphics, images, and data. These cryo inkjet labels also display excellent resistance against various lab chemicals, such as alcohols, acetone, xylene, as well as sanitizing wipes.



#CIJ-29



#CIJ-28

Sheet Format

Temperature range -196°C to +121°C (-321°F to +250°F)

Our sheet format Cryo-JetTAG labels can be printed with most standard desktop inkjet printers or written on using cryo markers and ballpoint pens. Perfect for identifying containers destined for prolonged storage in cryogenic conditions, including in liquid nitrogen Dewars (-196°C) and ultra-low temperature lab freezers (-80°C). Free online MS Word templates are available for quick and easy template set-up.



#CIJS-23



#CIJS-3



Cryo-JetTAG™ Inkjet Cryogenic Labels - Roll Format



#CIJ-29




#CIJ-28



#CIJ-29

Item #	Inches (W x H)	mm (W x H)	Perforation	Core	Labels/Roll
CIJ-277NPC1-1WH	0.625" x 1.2"	15.9 x 30.5	No	1"	1000
CIJ-277NPC2-1WH	0.625" x 1.2"	15.9 x 30.5	No	2"	1000
CIJ-277NPC3-1WH	0.625" x 1.2"	15.9 x 30.5	No	3"	1000
CIJ-279C1-1WH	1" x 1.5"	25.4 x 38.1	Yes	1"	1000
CIJ-279C2-1WH	1" x 1.5"	25.4 x 38.1	Yes	2"	1000
CIJ-279C3-3WH	1" x 1.5"	25.4 x 38.1	Yes	3"	3000
CIJ-183C1-2WH	1" x 0.5" + 0.375"	25.4 x 12.7 + 9.5	Yes	1"	2000
CIJ-183C2-2WH	1" x 0.5" + 0.375"	25.4 x 12.7 + 9.5	Yes	2"	2000
CIJ-183C3-2WH	1" x 0.5" + 0.375"	25.4 x 12.7 + 9.5	Yes	3"	2000
CIJ-29C1-1WH	1" x 1"	25.4 x 25.4	Yes	1"	1000
CIJ-29C2-1WH	1" x 1"	25.4 x 25.4	Yes	2"	1000
CIJ-29C3-1WH	1" x 1"	25.4 x 25.4	Yes	3"	1000
CIJ-538C1-1WH	1.125" x 1"	28.6 x 25.4	Yes	1"	1000
CIJ-538C2-1WH	1.125" x 1"	28.6 x 25.4	Yes	2"	1000
CIJ-538C3-2WH	1.125" x 1"	28.6 x 25.4	Yes	3"	2000
CIJ-10C1-1WH	1.5" x 0.75"	38.1 x 19.1	Yes	1"	1000
CIJ-10C2-1WH	1.5" x 0.75"	38.1 x 19.1	Yes	2"	1000
CIJ-10C3-1WH	1.5" x 0.75"	38.1 x 19.1	Yes	3"	1000
CIJ-28C1-2WH	2" x 1"	50.8 x 25.4	Yes	1"	2000
CIJ-28C2-2WH	2" x 1"	50.8 x 25.4	Yes	2"	2000
CIJ-28C3-2WH	2" x 1"	50.8 x 25.4	Yes	3"	2000
CIJ-40NPC1-1WH	2" x 1.125"	50.8 x 28.6	No	1"	1000
CIJ-40NPC2-1WH	2" x 1.125"	50.8 x 28.6	No	2"	1000
CIJ-40NPC3-2WH	2" x 1.125"	50.8 x 28.6	No	3"	2000

Cryo-JetTAG™ Inkjet Cryogenic Labels - Sheet Format



#CIJS-23

Item #	Inches (W x H)	mm (W x H)	Sheet Format	Labels/Sheet	Sheets/Pack
CIJS-23WH	1.28" x 0.5"	32.5 x 12.7	8.5" x 11"	126	16
CIJS-8WH	1.77" x 0.79"	45 x 20.1	8.5" x 11"	52	16
CIJS-3WH	2.63" x 1"	66.8 x 25.4	8.5" x 11"	30	16

Shop on www.Labtag.com/cryojettag