

-196°C



Cryo-LazrTAG™

Cryogenic Labels for Laser Printers

Labtag.com/cryo-lazrtag



Cryo-LazrTAG™

Cryogenic Labels for Laser Printers




Cryo-LazrTAG™ thermoplastic waterproof labels are specifically designed to adhere to vials, microtubes, microplates, glass tubes, polycarbonate cryo/freezer boxes, and other containers that will undergo long-term deep-freeze or cryogenic storage in liquid nitrogen, laboratory freezers, as well as transportation on dry ice. Resistant to multiple freeze-thaw cycles without fading, detaching, or wrinkling, these laser labels are provided in sheet format for printing with desktop laser printers, are wear-and-tear proof, will not jam printers, and resist UV and gamma irradiation.

Permanent | Hagaki Format (4" x 6")

Temperature range -196°C to +120°C (-321°F to +248°F)

Cryo-LazrTAG™ cryogenic plastic rectangle labels designed for 0.2 ml PCR tubes, 0.5 ml, 1.5 ml, and 2ml microcentrifuge tubes, and 1 ml to 2 ml cryo vials, that can withstand deep-freezing, PCR thermocycler temperatures, and multiple freeze-thaw cycles.




SKU#	Inches (W x H)	mm (W x H)	Available Colors
CLH-2	0.79" x 0.2"	20 x 5.1	WH SI BE YE LE GA MT PI LA BL OR BR RE VI GS LA
CLH-1	1.02" x 0.4"	26 x 10	WH SI BE YE LE GA MT PI LA BL OR BR RE VI GS LA

Permanent | US Letter (8.5" x 11")

Temperature range -196°C to +120°C (-321°F to +248°F)

Cryo-LazrTAG™ labels coated with a permanent adhesive for long-term cryostorage of cryo vials, tubes, microplates, and polycarbonate cryo/freezer boxes.




SKU#	Inches (W x H)	mm (W x H)	Colors Available
*HCL-10	0.35"	9	WH SI YE LE GA MT PI LA BL OR BR RE LA
*CL-48T1	0.433"	11	WH SI YE LE GA MT PI LA BL OR BR RE LA
CL-51T1	0.75"	19.1	WH LE
CL-16T1	0.75" x 0.75"	19.1 x 19.1	WH OR
CL-12T1	0.94" x 0.5"	23.9 x 12.7	WH SI YE LE GA MT PI LA BL OR RE GS LA
CL-32T1	0.94" x 0.7813"	23.9 x 19.8	WH GA
CL-71T1	1" x 1"	25.4 x 25.4	WH
CL-60T1	1.125" x 0.94" + 0.44"	28.6 x 23.9 + 11.2	WH
CL-59T1	1.2" x 0.63"	30.5 x 16	WH
CL-9T1	1.26" x 0.87"	32 x 22.1	WH SI BE YE LE GA MT PI LA BL OR RE LA
CL-37T1	1.26" x 1.1"	32 x 27.9	WH YE BL
CL-23T1	1.28" x 0.5"	32.5 x 12.7	WH SI BE YE LE GA MT PI LA BL OR RE LA
CL-47T1	1.3" x 1.3"	33 x 33	WH
CL-6T1	1.42" x 0.55"	36 x 14	WH SI BE YE LE GA MT PI LA BL OR RE LA
CL-17T1	1.5" x 0.25"	38.1 x 6.4	WH
CL-13T1	1.75" x 0.5"	44.5 x 12.7	WH SI BE YE LE GA MT PI LA BL OR BR RE LA
CL-8T1	1.77" x 0.79"	45 x 20	WH SI BE YE LE GA MT PI LA BL OR BR RE LA
CL-33T1	1.97" x 0.97"	50 x 24.6	WH GA RE
CL-63T1	2" x 1.125"	50.8 x 28.6	WH
CL-3T1	2.63" x 1"	66.7 x 25.4	WH YE GA PI BL OR RE GS LA
CL-29T1	2.83" x 1.57"	71.9 x 39.9	WH YE GA PI LA BL OR RE
CL-30T1	3" x 2.375"	76.2 x 60.3	WH YE
CL-7T1	4" x 2"	101.6 x 50.8	WH OR BR
CL-15T1	4" x 3.33"	101.6 x 84.6	WH YE GA PI LA BL OR

*HCL-10/CL-48: Deep-freeze labels for freezer storage (-80°C/-112°F)

Removable | US Letter (8.5" x 11")

Temperature range -196°C to +120°C (-321°F to +248°F)

Removable Cryo-LazrTAG™ are durable, glove-friendly thermoplastic laser labels for cryogenic applications, removable from all surfaces without leaving behind any residual adhesive.




SKU#	Inches (W x H)	mm (W x H)	Colors Available
RCL-12T1	0.94" x 0.5"	23.9 x 12.7	WH SI BE YE GA PI LA BL OR BR RE GS A
RCL-9T1	1.26" x 0.887"	32 x 22.1	WH SI BE YE GA PI LA BL OR RE GS A
RCL-37T1	1.26" x 1.1"	32 x 27.9	WH
RCL-23T1	1.28" x 0.5"	32.5 x 12.7	WH YE LE GA MT PI LA BL OR RE A
RCL-6T1	1.42" x 0.55"	36 x 14	WH SI BE YE LE GA MT PI LA BL OR PU RE VI GS A
RCL-17T1	1.5" x 0.25"	38.1 x 6.4	WH SI YE GA MT PI LA BL OR RE
RCL-13T1	1.75" x 0.5"	44.5 x 12.7	WH GA OR
RCL-76T1	1.75" x 2"	44.5 x 50.8	WH
RCL-8T1	1.77" x 0.79"	45 x 20	WH SI BE YE LE GA MT PI LA BL OR BR RE VI GS A
RCL-33T1	1.97" x 0.97"	50 x 24.6	WH GA RE
RCL-20T1	2" x 2"	50.8 x 50.8	WH
RCL-3T1	2.63" x 1"	66.8 x 25.4	WH SI BE YE LE GA MT PI LA BL OR RE GS A
RCL-11T1	2.64" x 0.277"	67.1 x 7	WH SI BE YE LE GA PI LA BL OR RE GS A
RCL-19T1	2.75" x 2.75"	69.9 x 69.9	WH YE BL
RCL-61T1	3.5" x 2"	88.9 x 50.8	WH
RCL-26T1	3.5" x 5"	88.9 x 127	WH
RCL-7T1	4" x 2"	101.6 x 50.8	WH SI BE YE GA PI LA BL OR BR RE GS
RCL-15T1	4" x 3.33"	101.6 x 84.6	WH YE GA PI LA BL OR RE GS

Wrap-Around Permanent | US Letter (8.5" x 11")

Temperature range -196°C to +150°C (-321°F to +302°F)

Self-laminating Cryo-LazrTAG™ are cryogenic laser labels for long-term storage in harsh ultra-low temperature environments. The clear laminating part of the labels protects the printed information from abrasion, and chemical exposure. These versatile labels can also withstand high-temperatures and autoclaving.




#DFSL-62

SKU#	Inches (W x H + Wrap)	mm (W x H + Wrap)	Labels/Pack	Colors Available
DFSL-13WH	1" x 0.5" + 0.75"	25.4 x 12.7 + 19.1	3,200	WH
DFSL-65	1" x 0.6" + 1.6"	25.4 x 15.2 + 40.6	1,920	WH, YE, MT, PI
DFSL-75WH	1" x 0.625" + 2"	25.4 x 15.9 + 50.8	1,280	WH
DFSL-77WH	2" x 0.625" + 2"	50.8 x 15.9 + 50.8	1,280	WH

Transparent Permanent | US Letter (8.5" x 11")

Temperature range -196°C to +150°C (-321°F to +302°F)

Transparent Cryo-LazrTAG™ cryogenic thermoplastic laser labels for long-term cryostorage of cryo vials, tubes, and boxes that allow the contents of labeled containers to remain clearly visible and give a "printed-on effect". Resistant to high-temperatures and autoclaving.



#DFLT-65


SKU#	Inches (W x H)	mm (W x H)	Labels/Pack
*DFLT-13	1.75" x 0.5"	44.5 x 12.7	1,280
DFLT-65	2.6" x 0.6"	66 x 15.2	768
DFLT-75	3" x 0.625"	76.2 x 15.9	512
DFLT-77	4" x 0.625"	101.6 x 15.9	512

*These transparent laser labels are also available with a color tab, in assorted colors.

Transparent Removable | US Letter (8.5" x 11")

Temperature range -196°C to +120°C (-321°F to +248°F)

Cryo-LazrTAG™ labels are waterproof and tear-proof, for the identification of flat surface containers in liquid nitrogen and laboratory freezers. With a glove-friendly adhesive, these transparent cryogenic labels can be removed without leaving any adhesive residue on the surface.



#TRCL-3

SKU#	Inches (W x H)	mm (W x H)	Labels/Pack
TRCL-12	0.94" x 0.5"	23.9 x 12.7	2,688
TRCL-23	1.28" x 0.5"	32.5 x 12.7	2,016
TRCL-13	1.75" x 0.5"	44.5 x 12.7	1,280
TRCL-3	2.63" x 1"	66.8 x 25.4	480

ABOUT US

GA International is a leading manufacturer of high-performance specialty labels, and a supplier of identification solutions used in research and clinical labs performing histological staining, cryo-preservation and sterilization. We understand the importance of proper patient and sample identification.

Designed by scientists for scientists, our labels endure the harshest laboratory applications to protect your valuable information.

Xylene Labels • Cryogenic Labels • Autoclave Labels • Healthcare Labels

Tapes • Ribbons • Printers • Scanners • Software • Automation

Printing Services • Data Printing • Industry Solutions • Custom Manufacturing

ISO 9001:2015 Certified



USA / CANADA
+1-800-518-0364

INTERNATIONAL
+1-450-973-9420

info@ga-international.com
www.labtag.com