



# Cryo-LazrTAG™

## Cryogenic Labels for Laser Printers



Cryo-LazrTAG™ thermoplastic waterproof labels are specifically designed to adhere to vials, microtubes, microplates, glass tubes, polycarbonate cryo/freezer boxes, and other containers that will undergo long-term deep-freeze or cryogenic storage in liquid nitrogen, laboratory freezers, as well as transportation on dry ice. Resistant to multiple freeze-thaw cycles without fading, detaching, or wrinkling, these laser labels are provided in sheet format for printing with desktop laser printers, are wear-and-tear proof, will not jam printers, and resist UV and gamma irradiation.



[Labtag.com/cryo-lazrtag](http://Labtag.com/cryo-lazrtag)

# Permanent | US Letter and Hagaki Format

Temperature range -196°C to +120°C (-321°F to +248°F)

Our permanent Cryo-LazrTAG™ cryogenic labels are designed to identify cryo vials, microplates, cryo/freezer boxes, as well as PCR and microcentrifuge tubes for extended storage in liquid nitrogen and lab freezers. Thermocycler-safe, they will also withstand multiple freeze-thaw cycles.

## Hagaki Format (4" x 6")

SKU#	Inches (W x H)	mm (W x H)	Colors Available
CLH-2	0.79" x 0.2"	20 x 5.1	(WH) SI BE YE LE GA MT PI LA BL OR BR RE VI GS A
CLH-1	1.02" x 0.4"	26 x 10	(WH) SI BE YE LE GA MT PI LA BL OR BR RE VI GS A

## US Letter Format (8.5" x 11")

SKU#	Inches (W x H)	mm (W x H)	Colors Available
*HCL-10	0.35"	9	(WH) SI YE LE GA MT PI LA BL OR BR RE VI A
CL-48T1	0.433"	11	(WH) SI YE LE GA MT PI LA BL OR BR RE VI A
CL-51T1	0.75"	19.1	(WH) LE
CL-16T1	0.75" x 0.75"	19.1 x 19.1	(WH) OR
CL-12T1	0.94" x 0.5"	23.9 x 12.7	(WH) SI YE LE GA MT PI LA BL OR RE GS A
CL-32T1	0.94" x 0.7813"	23.9 x 19.8	(WH) GA
CL-71T1	1" x 1"	25.4 x 25.4	(WH)
CL-60T1	1.125" x 0.94" + 0.44"	28.6 x 23.9 + 11.2	(WH)
CL-59T1	1.2" x 0.63"	30.5 x 16	(WH)
CL-9T1	1.26" x 0.87"	32 x 22.1	(WH) SI BE YE LE GA MT PI LA BL OR RE A
CL-37T1	1.26" x 1.1"	32 x 27.9	(WH) YE BL
CL-23T1	1.28" x 0.5"	32.5 x 12.7	(WH) SI BE YE LE GA MT PI LA BL OR RE A
CL-47T1	1.3" x 1.3"	33 x 33	(WH)
CL-6T1	1.42" x 0.55"	36 x 14	(WH) SI BE YE LE GA MT PI LA BL OR RE A
CL-17T1	1.5" x 0.25"	38.1 x 6.4	(WH)
CL-13T1	1.75" x 0.5"	44.5 x 12.7	(WH) SI BE YE LE GA MT PI LA BL OR BR RE A
CL-8T1	1.77" x 0.79"	45 x 20	(WH) SI BE YE LE GA MT PI LA BL OR BR RE A
CL-33T1	1.97" x 0.97"	50 x 24.6	(WH) GA RE
CL-63T1	2" x 1.125"	50.8 x 28.6	(WH)
CL-3T1	2.63" x 1"	66.7 x 25.4	(WH) YE GA PI BL OR RE GS A
CL-29T1	2.83" x 1.57"	71.9 x 39.9	(WH) YE GA PI LA BL OR RE
CL-30T1	3" x 2.375"	76.2 x 60.3	(WH) YE
CL-7T1	4" x 2"	101.6 x 50.8	(WH) OR BR
CL-15T1	4" x 3.33"	101.6 x 84.6	(WH) YE GA PI LA BL OR

\*HCL-10/CL-48: Deep-freeze labels for freezer storage (-80°C/-112°F)




Available in A4 sheet format, with more sizes, colors and configurations on [Labtag.com/cryo-lazrtag](http://Labtag.com/cryo-lazrtag)

# Removable | US Letter (8.5" x 11")

Temperature range -196°C to +120°C (-321°F to +248°F)

Removable Cryo-LazrTAG™ are durable, glove-friendly thermoplastic laser labels for cryogenic applications, removable from all surfaces without leaving behind any residual adhesive.



SKU#	Inches (W x H)	mm (W x H)	Colors Available
RCL-12T1	0.94" x 0.5"	23.9 x 12.7	(WH) SI BE YE GA PI LA BL OR BR RE GS A
RCL-9T1	1.26" x 0.887"	32 x 22.1	(WH) SI BE YE GA PI LA BL OR RE GS A
RCL-37T1	1.26" x 1.1"	32 x 27.9	(WH)
RCL-23T1	1.28" x 0.5"	32.5 x 12.7	(WH) YE LE GA MT PI LA BL OR RE A
RCL-6T1	1.42" x 0.55"	36 x 14	(WH) SI BE YE LE GA MT PI LA BL OR PU RE VI GS A
RCL-17T1	1.5" x 0.25"	38.1 x 6.4	(WH) SI YE GA MT PI LA BL OR RE
RCL-13T1	1.75" x 0.5"	44.5 x 12.7	(WH) GA OR
RCL-76T1	1.75" x 2"	44.5 x 50.8	(WH)
RCL-8T1	1.77" x 0.79"	45 x 20	(WH) SI BE YE LE GA MT PI LA BL OR BR RE VI GS A
RCL-33T1	1.97" x 0.97"	50 x 24.6	(WH) GA RE
RCL-20T1	2" x 2"	50.8 x 50.8	(WH)
RCL-3T1	2.63" x 1"	66.8 x 25.4	(WH) SI BE YE LE GA MT PI LA BL OR RE GS A
RCL-11T1	2.64" x 0.277"	67.1 x 7	(WH) SI BE YE LE GA PI LA BL OR RE GS A
RCL-19T1	2.75" x 2.75"	69.9 x 69.9	(WH) YE BL
RCL-61T1	3.5" x 2"	88.9 x 50.8	(WH)
RCL-26T1	3.5" x 5"	88.9 x 127	(WH)
RCL-7T1	4" x 2"	101.6 x 50.8	(WH) SI BE YE GA PI LA BL OR BR RE GS
RCL-15T1	4" x 3.33"	101.6 x 84.6	(WH) YE GA PI LA BL OR RE GS



Available in A4 sheet format, with more sizes, colors and configurations on [Labtag.com/cryo-lazrtag](http://Labtag.com/cryo-lazrtag)

## Wrap-Around Permanent | US Letter (8.5" x 11")

Temperature range -196°C to +150°C (-321°F to +302°F)

Self-laminating Cryo-LazrTAG™ are cryogenic laser labels for long-term storage in harsh ultra-low temperature environments. The clear laminating part of the labels protects the printed information from abrasion, and chemical exposure. These versatile labels can also withstand high-temperatures and autoclaving.



SKU#	Inches (W x H + Wrap)	mm (W x H + Wrap)	Labels/Pack	Colors Available
DFSL-13WH	1" x 0.5" + 0.75"	25.4 x 12.7 + 19.1	3,200	WH
DFSL-65	1" x 0.6" + 1.6"	25.4 x 15.2 + 40.6	1,920	WH YE MT PI
DFSL-75WH	1" x 0.625" + 2"	25.4 x 15.9 + 50.8	1,280	WH
DFSL-77WH	2" x 0.625" + 2"	50.8 x 15.9 + 50.8	1,280	WH

## Transparent Permanent | US Letter (8.5" x 11")

Temperature range -196°C to +150°C (-321°F to +302°F)

Transparent Cryo-LazrTAG™ cryogenic thermoplastic laser labels for long-term cryostorage of cryo vials, tubes, and boxes that allow the contents of labeled containers to remain clearly visible and give a "printed-on effect". Resistant to high-temperatures and autoclaving.



SKU#	Inches (W x H)	mm (W x H)	Labels/Pack
*DFLT-13	1.75" x 0.5"	44.5 x 12.7	1,280
DFLT-65	2.6" x 0.6"	66 x 15.2	768
DFLT-75	3" x 0.625"	76.2 x 15.9	512
DFLT-77	4" x 0.625"	101.6 x 15.9	512

\*These transparent laser labels are also available with a color tab, in assorted colors.

## Transparent Removable | US Letter (8.5" x 11")

Temperature range -196°C to +120°C (-321°F to +248°F)

Cryo-LazrTAG™ labels are waterproof and tear-proof, for the identification of flat surface containers in liquid nitrogen and laboratory freezers. With a glove-friendly adhesive, these transparent cryogenic labels can be removed without leaving any adhesive residue on the surface.



SKU#	Inches (W x H)	mm (W x H)	Labels/Pack
TRCL-12	0.94" x 0.5"	23.9 x 12.7	2,688
TRCL-23	1.28" x 0.5"	32.5 x 12.7	2,016
TRCL-13	1.75" x 0.5"	44.5 x 12.7	1,280
TRCL-3	2.63" x 1"	66.8 x 25.4	480