

Cryo-WrapTAG™


Wrap-Around Cryogenic Thermal-Transfer Labels



- Withstand long-term storage under cryogenic conditions, including vapor and liquid phase liquid nitrogen (-196°C/-321°F)
- Our HBTT-class wrap-around labels can also resist elevated temperatures, up to 100°C/212°F
- CATT-class wrap-around labels can resist temperatures up to 150°C/302°F, able to withstand sterilization protocols
- Self-laminating, the wrap-around clear laminate protects the printout from abrasion and scratching
- Able to withstand exposure to harsh chemicals and solvents, such as alcohol, bleach and detergents
- Strong permanent adhesive, well suited for small diameter cryo vials and microtubes

Recommended Ribbons: RR-class

Shop on www.Labtag.com/cryo-wraptag



Self-Laminating Area

Printable Area



Cryogenic

Temperature range -196°C to +100°C (-321°F to +212°F)

SKU#	Inches (W x H)	mm (W x H)	Core	Perforation	Labels/Roll	Colors Available
HBTT-337NPNOTC1-1WH	0.59" x 1.06" + 1.38"	15 x 26.9 + 35.1	1"	No	1,000	(WH)
HBTT-337NPNOTC3-1WH	0.59" x 1.06" + 1.38"	15 x 26.9 + 35.1	3"	No	1,000	(WH)
HBTT-338NPNOTC1-1WH	0.59" x 1.06" + 0.59"	15 x 26.9 + 15	1"	No	1,000	(WH)
HBTT-338NPNOTC3-2WH	0.59" x 1.06" + 0.59"	15 x 26.9 + 15	3"	No	2,000	(WH)
HBTT-344NOTC1-1WH	0.866" x 0.866" + 0.866"	22 x 22 + 22	1"	Yes	1,000	(WH)
HBTT-344NOTC3-2WH	0.866" x 0.866" + 0.866"	22 x 22 + 22	3"	Yes	2,000	(WH)
HBTT-311NOTC1-2WH	1" x 0.625" + 1"	25.4 x 15.9 + 25.4	1"	Yes	2,000	(WH)
HBTT-311NOTC3-4WH	1" x 0.625" + 1"	25.4 x 15.9 + 25.4	3"	Yes	4,000	(WH)
HBTT-312NOTC1-1WH	1" x 0.625" + 1.375"	25.4 x 15.9 + 34.9	1"	Yes	1,000	(WH)
HBTT-312NOTC3-3WH	1" x 0.625" + 1.375"	25.4 x 15.9 + 34.9	3"	Yes	3,000	(WH)
HBTT-322NOTC1-2WH	1" x 0.75" + 1.75"	25.4 x 19.1 + 44.5	1"	Yes	2,000	(WH)
HBTT-322NOTC3-4WH	1" x 0.75" + 1.75"	25.4 x 19.1 + 44.5	3"	Yes	4,000	(WH)
HBTT-317NOTC1-1WH	1" x 1" + 1.72"	25.4 x 25.4 + 43.7	1"	Yes	1,000	(WH)
HBTT-317NOTC3-2WH	1" x 1" + 1.72"	25.4 x 25.4 + 43.7	3"	Yes	2,000	(WH)
HBTT-309NOTC1-1	1" x 1.18" + 1.69"	25.4 x 30 + 42.9	1"	Yes	1,000	(WH) (YE) (GA) (PI)
HBTT-309NOTC3-1	1" x 1.18" + 1.69"	25.4 x 30 + 42.9	3"	Yes	1,000	(WH) (YE) (GA) (PI)
HBTT-339NPNOTC1-1	1.06" x 0.59" + 0.59"	26.9 x 15 + 15	1"	No	1,000	(WH) (YE) (GA) (OR) (BL) (LA) (SI)
HBTT-339NPNOTC3-2	1.06" x 0.59" + 0.59"	26.9 x 15 + 15	3"	No	2,000	(WH) (YE) (GA) (OR) (BL) (LA) (SI)
HBTT-303NOTC1-1	1.125" x 1" + 1.595"	28.6 x 25.4 + 40.5	1"	Yes	1,000	(WH) (YE) (GA) (PI) (OR) (RE) (BL) (LA) (SI)
HBTT-303NOTC3-1	1.125" x 1" + 1.595"	28.6 x 25.4 + 40.5	3"	Yes	1,000	(WH) (YE) (GA) (PI) (OR) (RE) (BL) (LA) (SI)

Cryogenic & Autoclave-Resistant

Temperature range -196°C to +150°C (-321°F to +302°F)

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll	Colors Available
CATT-318SBNOTC1-1WH	1" x 0.25" + 0.75"	25.4 x 12.7 + 19.1	1"	1,000	(WH)
CATT-318SBNOTC3-2WH	1" x 0.25" + 0.75"	25.4 x 12.7 + 19.1	3"	2,000	(WH)
CATT-312NOTC1-1WH	1" x 0.625" + 1.375"	25.4 x 15.9 + 35	1"	1,000	(WH)
CATT-312NOTC3-2WH	1" x 0.625" + 1.375"	25.4 x 15.9 + 35	3"	2,000	(WH)
CATT-311NOTC1-2WH	1" x 0.625" + 1"	25.4 x 15.9 + 25.4	1"	2,000	(WH)
CATT-311NOTC3-4WH	1" x 0.625" + 1"	25.4 x 15.9 + 25.4	3"	4,000	(WH)
CATT-309NOTC1-1	1" x 1.18" + 1.69"	25.4 x 30 + 42.9	1"	1,000	(WH) (YE) (GA) (PI)
CATT-309NOTC3-1	1" x 1.18" + 1.69"	25.4 x 30 + 42.9	3"	1,000	(WH) (YE) (GA) (PI)
CATT-319NOTC1-1	1.18" x 1" + 0.75"	30 x 25.4 + 19.1	1"	1,000	(WH) (YE) (GA) (PI)
CATT-319NOTC3-2	1.18" x 1" + 0.75"	30 x 25.4 + 19.1	3"	2,000	(WH) (YE) (GA) (PI)
CATT-303NOTC1-1WH	1.125" x 1" + 1.595"	28.6 x 25.4 + 40.5	1"	1,000	(WH)
CATT-303NOTC3-1WH	1.125" x 1" + 1.595"	28.6 x 25.4 + 40.5	3"	1,000	(WH)
CATT-308NOTC1-1WH	1.34" x 1" + 1.38"	34 x 25.4 + 35	1"	1,000	(WH)
CATT-308NOTC3-1WH	1.34" x 1" + 1.38"	34 x 25.4 + 35	3"	1,000	(WH)