

Cryo-WrapTAG™

Wrap-Around Thermal-Transfer Labels



These self-laminating, wrap-around cryogenic labels have a clear laminate that is applied over the printed portion of the label, protecting it against harsh laboratory conditions, such as abrasion, fading, exposure to harsh chemicals, and extreme temperatures. They also withstand long-term storage under cryogenic conditions, including in liquid phase and vapor phase liquid nitrogen (-196°C/-321°F).



Cryogenic Wrap-Around Labels

Temperature range -196°C to +100°C (-321°F to +212°F)

Our HBTT-class labels are designed for the secure identification of small diameter cryo vials and microtubes.



SKU#	Inches (W x H)	mm (W x H)	Available Cores	Perforation	Labels/Roll	Colors Available
HBTT-338NPNOT	0.59" x 1.06" + 0.59"	15 x 26.9 + 15	1", 3"	No	1,000, 2,000	(WH)
HBTT-344NOT	0.866" x 0.866" + 0.866"	22 x 22 + 22	1", 3"	Yes	1,000, 2,000	(WH)
HBTT-311NOT	1" x 0.625" + 1"	25.4 x 15.9 + 25.4	1", 3"	Yes	2,000, 4,000	(WH)
HBTT-312NOT	1" x 0.625" + 1.375"	25.4 x 15.9 + 34.9	1", 3"	Yes	1,000, 3,000	(WH)
HBTT-317NOT	1" x 1" + 1.72"	25.4 x 25.4 + 43.7	1", 3"	Yes	1,000, 2,000	(WH)
HBTT-309NOT	1" x 1.18" + 1.69"	25.4 x 30 + 42.9	1", 3"	Yes	1,000	(WH) YE GA PI
HBTT-339NPNOT	1.06" x 0.59" + 0.59"	26.9 x 15 + 15	1", 3"	No	1,000, 2,000	(WH) YE GA OR BL LA SI
HBTT-303NOT	1.125" x 1" + 1.595"	28.6 x 25.4 + 40.5	1", 3"	Yes	1,000	(WH) YE GA PI OR RE BL LA SI



Cryogenic & Autoclave-Resistant Wrap-Around Labels

Temperature range -196°C to +150°C (-321°F to +302°F)

Our CATT-class wrap-around labels are also able to withstand elevated temperatures, utilized in various sterilization protocols.

SKU#	Inches (W x H)	mm (W x H)	Available Cores	Labels/Roll	Colors Available
CATT-318SBNOT	1" x 0.25" + 0.75"	25.4 x 12.7 + 19.1	1", 3"	1,000, 2,000	(WH)
CATT-312NOT	1" x 0.625" + 1.375"	25.4 x 15.9 + 35	1", 3"	1,000, 2,000	(WH)
CATT-311NOT	1" x 0.625" + 1"	25.4 x 15.9 + 25.4	1", 3"	2,000, 4,000	(WH)
CATT-309NOT	1" x 1.18" + 1.69"	25.4 x 30 + 42.9	1", 3"	1,000	(WH) YE GA PI
CATT-303NOT	1.125" x 1" + 1.595"	28.6 x 25.4 + 40.5	1", 3"	1,000	(WH)
CATT-308NOT	1.34" x 1" + 1.38"	34 x 25.4 + 35	1", 3"	1,000	(WH)




RFID Cryogenic Wrap-Around Labels (Patented)

Temperature Range -196°C to +100°C (-321°F to +212°F)

These labels combine RFID technology and the reliable wrap-around feature to ensure your cryo vials are securely identified with visual and electronic information.

SKU#	Inches (W x H)	mm (W x H)	Available Cores	Labels/Roll	Colors Available
RF3-HBTT-308NP	1.34" x 1" + 1.38"	34 x 25.4 + 35	1"	1.000	(WH)



Shop on www.Labtag.com/cryo-wraptag