

CryoSTUCK[®]

Cryo Labels for Frozen Surfaces

Temperature range from -196°C to +100°C (-321°F to + 212°F)

- Adhere to previously frozen vials, as low as -80°C (no thawing required)
- Resist extreme cryogenic temperatures (-196°C/-321°F)
- Ideal for over-labeling (blackout and clear options)
- Available for thermal-transfer and laser printers
- Perfect for biobanks, biorepositories, and cryogenic storage facilities
- Also available in tapes

Standard - Sheet Format



Provided in sheets for desktop laser printers, these permanent cryo labels will adhere to already frozen vials and containers. Free online templates make for easy setup and printing.

Standard - Roll Format



Cryogenic labels that adhere to already frozen surfaces, without the need to thaw them first. Available in clear and blackout formats, for printing with thermal-transfer barcode printers.

Blackout - Roll Format



The unique opaque design of these blackout CryoSTUCK labels allows them to cover-up existing labels. Perfect for relabeling frozen vials when removal of the current label is impossible.

Transparent - Roll Format



Our Clear CryoSTUCK labels provide a way to over-label cryo vials and tubes. This protects existing labels from harsh environments, without concealing their printed information.



Standard - Sheet Format (US Letter)

SKU#	W x H (Inches)	Labels/Sheet	Labels/Pack	Sheets/pack
LCS-9WH	1.26" x 0.87"	72	1,152	16
LCS-12WH	0.94" x 0.5"	168	2,688	16

Blackout - Roll Format

SKU#	W x H (Inches)	Core	Labels/Roll
AEA-5NOTC1-1WH	0.875" x 1.125" + 0.75"	1"	1,000
AEA-5NOTC3-1WH	0.875" x 1.125" + 0.75"	3"	1,000
AEA-7C1-1WH	1.57" x 0.75" + 0.93"	1"	1,000
AEA-7C3-2WH	1.57" x 0.75" + 0.93"	3"	2,000
L2FC-5C1-1WH	1.625" x 2"	1"	1,000
L2FC-5C3-1WH	1.625" x 2"	3"	1,000
AEA-3C1-1WH	2" x 1" + 0.5"	1"	1,000
AEA-3C3-1WH	2" x 1" + 0.5"	3"	1,000
L2FC-11C1-1WH	3" x 1.25"	1"	1,000
L2FC-11C3-1WH	3" x 1.25"	3"	1,000

Standard - Roll Format

SKU#	W x H (Inches)	Core	Labels/Roll
L2FS-5C1-1WH	0.875" x 2"	1"	1,000
L2FS-5C3-1WH	0.875" x 2"	3"	1,000
L2FS-1C1-1WH	1" x 0.5"	1"	1,000
L2FS-1C3-1WH	1" x 0.5"	3"	1,000
FIX-312NOTC1-1WH	1" x 0.625" + 1.375" wrap	1"	1,000
FIX-312NOTC3-3WH	1" x 0.625" + 1.375" wrap	3"	3,000
L2FS-2C1-2WH	1" x 1"	1"	2,000
L2FS-2C3-4WH	1" x 1"	3"	4,000
L2FS-4C1-2WH	1.25" x 0.875"	1"	2,000
L2FS-4C3-2WH	1.25" x 0.875"	3"	2,000
FIX-329NOTC1-1WH	1.5" x 1" + 1.5" wrap	1"	1,000
FIX-329NOTC3-2WH	1.5" x 1" + 1.5" wrap	3"	2,000
L2FS-3C1-1WH	3" x 0.55"	1"	1,000
L2FS-3C3-1WH	3" x 0.55"	3"	1,000

Transparent - Roll Format

SKU#	W x H (Inches)	Core	Labels/Roll
AHA-268NOTC1-2	1.75" x 0.75"	1"	2,000
AHA-268NOTC3-4	1.75" x 0.75"	3"	4,000
AHA-224NOTC1-2	1.75" x 1"	1"	2,000
AHA-224NOTC3-2	1.75" x 1"	3"	2,000
AHA-501NOTC1-0.25	3" x 4.75"	1"	250
AHA-501NOTC3-0.5	3" x 4.75"	3"	500



 [Shop on Labtag.com/cryostuck](https://www.Labtag.com/cryostuck)