

# DYMO-Compatible CryoSTUCK®

Patented

## Cryogenic Labels for Frozen Vials

Temperature range -196°C to +70°C (-321°F to +158°F)



Adhere as low as -80°C (-112°F)

These DYMO-compatible, solvent-resistant cryogenic labels will rapidly adhere to already frozen surfaces, such as cryo vials and tubes. They can be applied under deep-freeze conditions (as low as -80°C), eliminating the need to thaw samples beforehand, and can then be placed immediately back into storage, inside freezers (-80°C) and liquid nitrogen tanks (-196°C), with no cure time required. They display excellent resistance to laboratory chemicals like alcohols (100% Ethanol, immersion: up to 15 minutes), xylene (immersion: up to 30 minutes), detergents, and disinfecting solutions.



| SKU#       | Inches (W x H)               | mm (W x H)        | Labels/Roll |
|------------|------------------------------|-------------------|-------------|
| EDCS-040WH | 0.5" x 1"                    | 12.7 x 25.4       | 1,000       |
| EDCS-070WH | 1" x 1"                      | 25.4 x 25.4       | 750         |
| EDCS-081WH | 1.04" x 0.5" + 0.375" circle | 26.4 x 12.7 + 9.5 | 500 sets    |