

DiTi™ 203iD

The New Direct Thermal Printer by GA International

NEW!

Small, Lightweight, & Affordable

The new direct thermal printer from GA International, **DiTi 203iD**, can print our growing line of direct thermal labels designed for various applications throughout the lab, without the need for additional ink sources. Small and lightweight, the **DiTi 203iD** printer can be easily moved and set up anywhere in the lab, with an easy loading design for quick roll replacement.

- Accepts label rolls up to 3.54" wide
- Free barcoding software for designing labels
- 203 dpi print resolution
- USB



Technical Specifications

- **Print Method:** Direct Thermal
- **Print Resolution:** 203 dpi
- **Interfaces:** USB, Serial
- **Print Speed:** up to 4" (102mm/s)
- **Max. Print Width:** 3.35" (85.1mm)
- **Max. Print Length:** 11.8" (300mm)
- **Label Width:** 1" to 3.54" (25.4 to 90mm)
- **Label Height:** 0.59"-11.8" (15-300mm)
- **Media Types:** Continuous, Die-cut/gap Labels (Black Marked, Notch), Fanfold
- **Media Roll Core Size:** 1"
- **Media Maximum Roll OD:** 3.54" (90mm)

What's in the Box

- DiTi™ 203iD printer
- Power cord*
- Power adapter
- USB cable
- Flash drive with software and printer driver installation instructions



*The product includes a North American, European, or UK power cord and adapter. Other country-specific power cords must be purchased locally and separately.

Shop on [Labtag.com/diti-203id-printer](https://www.Labtag.com/diti-203id-printer)