



Complete Labeling Solutions for IVF & ART Facilities

Cryogenic Labels for Vitrification Devices

- Cryo-StrawTAG™ labels are compatible with all standard and high-security straws, as well as vitrification devices such as Cryotop®, and CryoLock®
- These labels are designed to withstand long-term cryo preservation inside liquid nitrogen (-196°C/-321°F), and transport on dry ice
- The cryo IVF straw labels withstand freezing and the vitrification of sperm, oocytes and embryos for long-term storage



SKU#	Inches (W x H)	mm (W x H)	Printable Area (Inches)	Printable Area (mm)	Core	Labels/Roll
CST-29NPNOT0.25C1-1WH	1" x 1"	25.4 x 25.4	1" x 0.25"	25.4 x 6.4	1"	1,000
CST-29NPNOT0.5C1-1WH	1" x 1"	25.4 x 25.4	1" x 0.375"	25.4 x 9.5	1"	1,000
CST-28NPNOT0.25C1-1WH	2" x 1"	50.8 x 25.4	2" x 0.25"	50.8 x 6.4	1"	1,000
CST-28NPNOT0.5C1-1WH	2" x 1"	50.8 x 25.4	2" x 0.375"	50.8 x 9.5	1"	1,000
CST-79NPNOT0.25C1-1WH	1.375" x 0.875"	35 x 22.2	1.375" x 0.25"	35 x 6.4	1"	1,000
CST-79NPNOT0.5C1-1WH	1.375" x 0.875"	35 x 22.2	1.375" x 0.375"	35 x 9.6	1"	1,000

Labeling Tubes & Vials

- Large catalog of cryogenic labels, designed for various applications and environments
- Thermal-transfer printable, available in multiple sizes, colors, and configurations
- Perfect for identifying cryo vials, microtubes, cryo/freezer boxes and other cryogenic containers

Standard Cryogenic Labels

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll
JTTA-7C1-1WH	1" x 0.5"	25.4 x 12.7	1"	1,000
JTTA-29C1-1WH	1" x 1"	25.4 x 25.4	1"	1,000
JTTA-9C1-1WH	1.25" x 0.875"	31.8 x 22.2	1"	1,000



Self-Laminating Cryo Labels for Added Protection

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll
HBTT-312NOTC1-1WH	1" x 0.625" + 1.375" wrap	25.4 x 15.9 + 34.9 wrap	1"	1,000
HBTT-317NOTC1-1WH	1" x 1" + 1.72" wrap	25.4 x 25.4 + 43.7 wrap	1"	1,000



Cryogenic Labels for Frozen Surfaces

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll
UC-98C1-1WH	0.875" x 0.875"	22.2 x 22.2	1"	1,000
UC-7C1-1WH	1" x 0.5"	25.4 x 12.7	1"	1,000
UC-52NOTC1-2WH	0.94" x 0.5" + 0.437"	23.9 x 12.7 + 11.1	1"	2,000
UC-29C1-1WH	1" x 1"	25.4 x 25.4	1"	1,000



Labeling IVF Culture Plates & Dishes

Labels for the Side of Plates

- Narrow labels made of thermoplastic film are perfect for labeling the sides of plates
- These labels are made to resist various laboratory environments, including incubators
- Plates with barcodes can be easily scanned with minimal manipulation needed

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll
FJT-209C1-4WH	1" x 0.25"	25.4 x 6.4	1"	4,000
FJT-169C1-4WH	1.25" x 0.25"	31.8 x 6.4	1"	4,000



Labels for the Bottom of Plates

- Transparent labels, with mirror-image printed text and barcode, are easily read and scanned through the plate when applied to the underside
- Quickly read and scan while your culture plate is still on the work area
- These labels are made to resist various laboratory environments, including incubators

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll
HBCL-254NOTC1-2	0.875" x 0.433"	22.2 x 11	1"	2,000
HBCL-7NPNOTC1-2	1" x 0.5"	25.4 x 12.7	1"	2,000



Labels for Cryo Canes, Metal Racks & Canisters

- Designed specifically to stick to metal surfaces like metal racks, aluminum cryo canes, and steel canisters, to be stored in cryogenic conditions, such as liquid nitrogen (-196°C/-321°F)
- Provides bird's-eye-view of the contents of your containers, eliminating the need to take containers out of storage
- Compatible with all common thermal-transfer printer brands


SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll
AWA-236C1-4WH	0.31" x 0.5"	7.9 x 12.7	1"	4,000
AWA-28C1-1WH	2" x 1"	50.8 x 25.4	1"	1,000
AWA-36C1-1WH	3" x 1"	76.2 x 25.4	1"	1,000
AWA-122C1-0.3WH	4" x 4"	101.6 x 101.6	1"	300



Patient Record & Datasheet Labels

We provide a wide range of paper labels, for barcode thermal-transfer printers. These labels can be printed with a crisp, clear printout, to ensure the scannability of your barcodes. The general use labels are also available in colors, for office use, patient records and requisition forms.

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll
GPA-7C1-3WH	1" x 0.5"	25.4 x 12.7	1"	3,000
GPA-28C1-1WH	2" x 1"	50.8 x 25.4	1"	1,000
GPA-127C1-1WH	4" x 1"	102 x 25.4	1"	1,000



Printers, Scanners & Software

Equip your lab with a complete labeling solution. LabTAG offers printers, barcoding software, scanners, everything you'll need to securely and accurately identify your files, samples and equipment.

Ready-made printing kits are available, or you can build your own IVF printing solution with the help of our scientific experts.



SKU#	Labels	Ribbon	Printer Features	Software	Scanner
SYS-ZD621-31-1 without scanner	3,000 Cryo-StrawTAG™ labels of your choice	1 Smudge-Proof Resin Ribbon 244'	USB, Serial, Ethernet, Bluetooth Interfaces	BarTender Professional Version software (1 Printer Edition)	-
SYS-ZD621-31-2 +scanner	3,000 Cryo-StrawTAG™ labels of your choice	1 Smudge-Proof Resin Ribbon 244'	USB, Serial, Ethernet, Bluetooth Interfaces	BarTender Professional Version software (1 Printer Edition)	EVO-S61-ADU-K



ISO 9001:2015 Certified
ISO 22301:2012 Compliant

+1-450-973-9420
cx@ga-international.com
www.Labtag.com

USA
P.O. Box 2633,
Champlain, NY
12919, United States
1-800-518-0364

CANADA
4455 Rue Louis-B.-Mayer
Laval QC,
H7P 6B5, Canada
1-800-518-0364

NETHERLANDS
Runweg 28A,
5258 BN Berlicum
Netherlands
+31 (0) 73-7370185+

UK
10 John Street,
London WC1N 2EB
United Kingdom
44 (20) 3769 5683