

Identification Solutions for Metal Surfaces

MetaliTAG™ cryogenic labels for the identification of metal surfaces, including aluminum cryogenic vial canes, metal racks for cryo and freezer boxes, as well as aluminum cassettes and metal sleeves. Specifically designed for easy marking and identification inside Dewar-type liquid nitrogen tanks (-196°C/-320.8°F) and laboratory freezers (-80°C/-112°F).

These specialty labels are made from a cryogenic material that has been tested and approved for applying to stainless steel and aluminum with subsequent immersion and storage in liquid phase or vapor phase liquid nitrogen. The MetaliTAG cryo labels can be produced in different colors and configurations upon request.



Laser MetaliTAG™

Temperature range from -196°C to +120°C (-320.8°F to +243.8°F)

Laser-printable MetaliTAG™ sheet-format labels are uniquely designed to adhere securely to metal racks and be stored in cryogenic conditions. They can be easily printed using standard office and desktop laser printers, providing high-quality and long-lasting results.

SKU#	Inches (W x H)	mm (W x H)	Sheet Size	Labels/Pack
AWL-63WH	2" x 1.125"	50.8 x 28.6	US Letter	432



MetaliTAG™ - Transparent

Temperature range from -196°C to +100°C (-320.8°F to +212°F)

Transparent, thermal transfer-printable MetaliTAG cryogenic labels are ideal for labeling tools and equipment with unique designs and logos, providing a "printed-on" effect.

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll
AWC-28NOTC1-1	2" x 1"	50.8 x 25.4	1"	1,000
AWC-28NOTC3-3	2" x 1"	50.8 x 25.4	3"	3,000




MetaliTAPE™

Temperature range from -196°C to +100°C (-320.8°F to 212°F)

MetaliTAPE is writable and designed for use on metal surfaces for long-term cryogenic storage. This tape should be applied at room temperature and can be used to identify stainless-steel and aluminum racks, X-ray cassettes, film cans, and more.

SKU#	Inches x Feet	mm x m	Core	Colors
TWB-19C1-50	0.75" x 50'	19mm x 15.24m	1"	WH YE GA BL OR RE
TWB-19C3-100	0.75" x 50'	19mm x 15.24m	3"	WH YE GA BL OR RE




MetaliTAG™ - White & Colors

Temperature range from -196°C to +100°C (-320.8°F to +212°F)

MetaliTAG cryogenic labels have been tested for stainless steel and aluminum with subsequent immersion and storage inside liquid or vapor phase liquid nitrogen.

AWB-232C1-2WH	Inches (W x H)	mm (W x H)	Core	Labels/ Roll	Colors Available
AWB-236C1-4WH	0.31" x 0.5"	7.9 x 12.7	1"	4,000	(WH)
AWB-236C3-8WH	0.31" x 0.5"	7.9 x 12.7	3"	8,000	(WH)
AWB-4C1-1WH	1" x 0.75"	25.4 x 19.1	1"	1,000	(WH)
AWB-4C3-4WH	1" x 0.75"	25.4 x 19.1	3"	4,000	(WH)
AWB-29C1-2WH	1" x 1"	25.4 x 25.4	1"	2,000	(WH)
AWB-29C3-2WH	1" x 1"	25.4 x 25.4	3"	2,000	(WH)
AWB-94C1-0.5WH	1" x 2"	25.4 x 50.8	1"	500	(WH)
AWB-94C3-2WH	1" x 2"	25.4 x 50.8	3"	2,000	(WH)
AWB-48NPC1-1WH	1.5" x 2.5"	38.1 x 63.5	1"	1,000	(WH)
AWB-48NPC3-1WH	1.5" x 2.5"	38.1 x 63.5	3"	1,000	(WH)
AWB-232C1-2WH	1.625" x 0.625"	41.3 x 15.9	1"	2,000	(WH)
AWB-232C3-4WH	1.625" x 0.625"	41.3 x 15.9	3"	4,000	(WH)
AWB-38C1-2WH	1.75" x 0.5"	44.5 x 12.7	1"	2,000	(WH)
AWB-38C3-6WH	1.75" x 0.5"	44.5 x 12.7	3"	6,000	(WH)
AWB-24C1-2WH	1.75" x 1"	44.5 x 25.4	1"	2,000	(WH)
AWB-24C3-2WH	1.75" x 1"	44.5 x 25.4	3"	2,000	(WH)
AWB-28C1-1	2" x 1"	50.8 x 25.4	1"	1,000	(WH) LYE MT OR YE RE VI BL GA
AWB-28C3-1	2" x 1"	50.8 x 25.4	3"	1,000	(WH) LYE MT OR YE RE VI BL GA
AWB-36C1-1WH	3" x 1"	76.2 x 25.4	1"	1,000	(WH)
AWB-36C3-1WH	3" x 1"	76.2 x 25.4	3"	1,000	(WH)
AWB-56C1-1WH	3" x 2"	76.2 x 50.8	1"	1,000	(WH)
AWB-56C3-1WH	3" x 2"	76.2 x 50.8	3"	1,000	(WH)
AWB-122C1-0.3WH	4" x 4"	101.6 x 101.6	1"	300	(WH)
AWB-122C3-0.5WH	4" x 4"	101.6 x 101.6	3"	500	(WH)



Cryo-ID-Kit™

Temperature range from -196°C to +80°C (-320.8°F to +176°F)

For the identification of metal racks and other containers for long-term storage in liquid nitrogen and ultra-low temperature freezers. Each Cryo-ID-Kit includes 12 plastic Cryo-Cards™ with connection hooks and writable cryo labels.

SKU#	Inches (W x H)	mm (W x H)	Tags/Pack	Colors/Pack
CID-02WH	2.13" x 1.18"	54.1 x 30	12	(WH)
CIDB-02YE	2.13" x 1.18"	54.1 x 30	12	(YE)
CIDB-02A	2.13" x 1.18"	54.1 x 30	12 (4 per color)	(YE) (MT) (OR) (BL)
CIDB-02A2	2.13" x 1.18"	54.1 x 30	12 (4 per color)	(YE) (MT) (OR) (RE)
CIDB-02A3	2.13" x 1.18"	54.1 x 30	12 (2 per color)	(YE) (MT) (OR) (RE) (VI) (BL) (GA)

