


XyliTAG™


Xylene & Chemical Resistant Labels



- Label, adhesive and ink withstand immersion in xylene, toluene, EtOH, IPA, methanol, acetone, hexane, DMSO, MEK, ACN, fixatives, histological stains, buffers and other chemicals.
- Intended for the identification of microscope slides, tissue cassettes and other flat surfaces.
- Waterproof, weatherproof, wear- and tear-proof.
- Available in numerous sizes.
- Supplied in roll format on 1" & 3" cores for printing with thermal transfer printers.

Recommended Ribbons: TC-Class Xylene Resistant Resin Ribbons.


Shop on www.Labtag.com/xyliTAG



Discover our products at www.Labtag.com





Xylene Resistant Labels



#HTT-117

0.866" x 0.59" / 22 x 15 mm



Item #	Core	Labels Across	Labels/Roll
HTT-117C1-3WH	1"	1	3,000
HTT-117C3-3WH	3"	1	3,000
Available Colors: <input type="radio"/> WH			



#HTT-133

0.866" x 0.75" / 22 x 19.1 mm



Item #	Core	Labels Across	Labels/Roll
HTT-133C1-4WH	1"	4	4,000
HTT-133C3-8WH	3"	4	8,000
Available Colors: <input type="radio"/> WH			



#HTT-98

0.875" x 0.875" / 22.2 x 22.2 mm



Item #	Core	Labels Across	Labels/Roll
HTT-98C1-1WH	1"	1	1,000
HTT-98C3-1WH	3"	1	1,000
Available Colors: <input type="radio"/> WH			



#HTT-152

0.9" x 0.75" / 23 x 19.1 mm


Item #	Core	Labels Across	Labels/Roll
HTT-152C1-1WH	1"	1	1,000
HTT-152C3-1WH	3"	1	1,000
Available Colors: <input type="radio"/> WH			



#HTT-67

0.9375" x 0.9375" / 23.8 x 23.8 mm

Item #	Core	Labels Across	Labels/Roll
HTT-67C1-1WH	1"	1	1,000
HTT-67C3-1WH	3"	1	1,000
Available Colors: <input type="radio"/> WH			



To request your **FREE** samples,
contact us at 1-800-518-0364 / 1-450-973-9420
support@ga-international.com
www.Labtag.com