

MetalITAG™

Cryogenic Labels for Metal Racks

Temperature range from -196°C to +100°C (-320.8°F to +212°F)



- Specifically designed to adhere to metal surfaces (e.g. stainless steel)
- Resist extreme cryogenic temperatures (-196°C/-320.8°F), in liquid or vapor phase liquid nitrogen)
- Printable with all common thermal-transfer barcode printer brands
- Versatile design allows labels to withstand temperatures up to +100°C/+212°F

MetalITAG™ cryogenic labels for the identification of metal surfaces in liquid nitrogen Dewars, including aluminum cryogenic vial canes, metal racks for cryo and freezer boxes, as well as aluminum cassettes and metal sleeves. Specifically designed for easy marking and identification inside Dewar-type liquid nitrogen tanks (-196°C/-320.8°F) and laboratory freezers (-80°C/-112°F).

These specialty labels are made from a cryogenic material that has been tested and approved for applying to stainless steel and aluminum with subsequent immersion and storage inside liquid phase or vapor phase liquid nitrogen. The MetalITAG™ cryo labels can be printed on with all common brands of thermal-transfer barcode printers and come in different colors and configurations upon request.




Shop on www.Labtag.com/metalitag


MetalTAG™ Cryogenic Labels for Metal Racks




Item #	Inches (W x H)	mm (W x H)	Core	Labels/Roll
AWA-29C1-2WH	1" x 1"	25.4 x 25.4	1"	2,000
AWA-29C3-2WH	1" x 1"	25.4 x 25.4	3"	2,000
Available Colors: (WH)				




Item #	Inches (W x H)	mm (W x H)	Core	Labels/Roll
AWA-48NPC1-1WH	1.5" x 2.5"	38.1 x 63.5	1"	1,000
AWA-48NPC3-1WH	1.5" x 2.5"	38.1 x 63.5	3"	1,000
Available Colors: (WH)				



Item #	Inches (W x H)	mm (W x H)	Core	Labels/Roll
AWA-28C1-1	2" x 1"	50.8 x 25.4	1"	1,000
AWA-28C3-1	2" x 1"	50.8 x 25.4	3"	1,000
Available Colors: (WH) (BE) (OR) (BR) (RE) (VI) (LYE) (LMT) (LPI) (LBL)				



Item #	Inches (W x H)	mm (W x H)	Core	Labels/Roll
AWA-36C1-1WH	3" x 1"	76.2 x 25.4	1"	1,000
AWA-36C3-1WH	3" x 1"	76.2 x 25.4	3"	1,000
Available Colors: (WH)				



Cryoid-Kit™ Hanging Tags

For the identification of metal racks and other containers for long-term storage in liquid nitrogen and ultra-low temperature freezers. Each Cryoid-Kit™ includes 12 plastic Cryo-Cards™ with connection hooks and thermal transfer cryo labels.



Item #	Inches (W x H)	mm (W x H)	Available Colors	Qty/Pk
CID-02WH	2.13" x 1.18"	54 x 30	White	12
CID-02YE	2.13" x 1.18"	54 x 30	Yellow	12
CID-02A	2.13" x 1.18"	54 x 30	Assorted 3 Colors*	12
CID-01WH	3.39" x 2.13"	86 x 54	White	12

*More colors available on www.Labtag.com

