

NitroTAG[®]

Thermal-Transfer Labels for Cryogenic Storage

Temperature range -196°C to +110°C (-321°F to +230°F)



- Resist extreme cryogenic temperatures (-196°C/-321°F, in liquid or vapor phase liquid nitrogen)
- Strong permanent cryo adhesive allows labels to adhere to wet surfaces
- Resist spraying and swabbing with alcohol and surface sanitizing wipes

Labtag.com/nitrotag



Rectangles

Rectangle NitroTAG labels are available in a range of sizes, perfect for identifying small diameter cryo vials, microtubes, cardboard or polypropylene cryo/freezer boxes, cell culture plates, and plastic bags. These thermal-transfer cryogenic labels have a strong permanent adhesive that will remain firmly attached, even in low temperature environments.



SKU#	Inches (W x H)	mm (W x H)	Core	Labels Across	Labels/Roll	Colors
JTTA-164C1-10	0.25" x 1"	6.4 x 25.4	1"	10	10,000	WH BL YE LA
JTTA-164C3-10	0.25" x 1"	6.4 x 25.4	3"	10	10,000	WH BL YE LA
JTTA-51C1-2WH	0.75" x 0.4"	19.1 x 10.2	1"	1	2,000	WH
JTTA-51C3-6WH	0.75" x 0.4"	19.1 x 10.2	3"	1	6,000	WH
JTTA-98C1-1WH	0.875" x 0.875"	22.2 x 22.2	1"	1	1,000	WH
JTTA-98C3-3WH	0.875" x 0.875"	22.2 x 22.2	3"	1	3,000	WH
JTTA-1C1-1	1" x 0.4375"	25.4 x 11.1	1"	1	1,000	WH PI LE
JTTA-1C3-6	1" x 0.4375"	25.4 x 11.1	3"	1	6,000	WH PI LE
JTTA-7C1-1	1" x 0.5"	25.4 x 12.7	1"	1	1,000	WH BL GA LE PI YE RE OR LA SI
JTTA-7C3-5	1" x 0.5"	25.4 x 12.7	3"	1	5,000	WH BL GA LE PI YE RE OR LA SI
JTTA-4C1-1WH	1" x 0.75"	25.4 x 19.1	1"	1	1,000	WH
JTTA-4C3-1WH	1" x 0.75"	25.4 x 19.1	3"	1	1,000	WH
JTTA-29C1-1	1" x 1"	25.4 x 25.4	1"	1	1,000	WH BL GA PI YE RE
JTTA-29C3-3	1" x 1"	25.4 x 25.4	3"	1	3,000	WH BL GA PI YE RE
JTTA-174C1-6WH	1" x 1"	25.4 x 25.4	1"	3	6,000	WH
JTTA-174C3-10WH	1" x 1"	25.4 x 25.4	3"	3	10,000	WH
JTTA-118C1-1WH	1" x 1.25"	25.4 x 31.8	1"	1	1,000	WH
JTTA-118C3-4WH	1" x 1.25"	25.4 x 31.8	3"	1	4,000	WH
JTTA-18C1-2WH	1.25" x 0.5"	31.8 x 12.7	1"	1	2,000	WH
JTTA-18C3-2WH	1.25" x 0.5"	31.8 x 12.7	3"	1	2,000	WH
JTTA-9C1-1	1.25" x 0.875"	31.8 x 22.2	1"	1	1,000	WH GA YE MT OR RE
JTTA-9C3-1	1.25" x 0.875"	31.8 x 22.2	3"	1	1,000	WH GA YE MT OR RE
JTTA-97C1-1WH	1.28" x 0.625"	32.5 x 15.9	1"	1	1,000	WH
JTTA-97C3-2WH	1.28" x 0.625"	32.5 x 15.9	3"	1	2,000	WH
JTTA-66C1-1	1.5" x 1.5"	38.1 x 38.1	1"	1	1,000	WH RE BGA
JTTA-66C3-2	1.5" x 1.5"	38.1 x 38.1	3"	1	2,000	WH RE BGA
JTTA-28C1-2	2" x 1"	50.8 x 25.4	1"	1	2,000	WH BL GA PI YE RE OR LA BR BL SI
JTTA-28C3-2	2" x 1"	50.8 x 25.4	3"	1	2,000	WH BL GA PI YE RE OR LA BR BL SI
JTTA-36C1-1	3" x 1"	76.2 x 25.4	1"	1	1,000	WH BL GA PI YE RE OR LA BR SI
JTTA-36C3-1	3" x 1"	76.2 x 25.4	3"	1	1,000	WH BL GA PI YE RE OR LA BR SI
JTTA-56C1-1	3" x 2"	76.2 x 50.8	1"	1	1,000	WH GA RE OR
JTTA-56C3-1	3" x 2"	76.2 x 50.8	3"	1	1,000	WH GA RE OR
JTTA-134C1-0.3WH	3" x 4"	76.2 x 101.6	1"	1	300	WH
JTTA-134C3-0.5WH	3" x 4"	76.2 x 101.6	3"	1	500	WH
JTTA-62C1-0.5	4" x 2"	101.6 x 50.8	1"	1	500	WH BL GA YE RE OR
JTTA-62C3-1	4" x 2"	101.6 x 50.8	3"	1	1,000	WH BL GA YE RE OR

Circles

NitroTAG dot labels are designed for the identification of vial caps and tube tops. Available in assorted colors, they are ideal for color coding sample containers. These cryo labels are writable but can also be printed with text and 2D barcodes for enhanced tracking using most common thermal-transfer printers.

SKU#	Inches (Diameter)	mm (Diameter)	Core	Labels Across	Labels/Roll	Colors
JTTA-59NPNOTC1-2	0.25"	6.4	1"	2	2,000	WH YE GA LE BL PI
JTTA-59NPNOTC3-2	0.25"	6.4	3"	2	2,000	WH YE GA LE BL PI
JTTA-61NPNOTC1-4	0.354"	9	1"	2	4,000	WH YE OR LA BL SI
JTTA-61NPNOTC3-4	0.354"	9	3"	2	4,000	WH YE OR LA BL SI
JTTA-22NOTC1-2	0.438"	11.1	1"	2	2,000	WH YE GA MT BL PI OR RE
JTTA-22NOTC3-4	0.438"	11.1	3"	2	4,000	WH YE GA MT BL PI OR RE
JTTA-500NOTC1-4	0.5"	12.7	1"	4	4,000	A Multicolor
JTTA-500NOTC3-8	0.5"	12.7	3"	4	8,000	A Multicolor
JTTA-293NOTC1-3	0.5"	12.7	1"	3	3,000	A Multicolor
JTTA-293NOTC3-6	0.5"	12.7	3"	3	6,000	A Multicolor
JTTA-53NOTC1-2	0.5"	12.7	1"	2	2,000	WH YE GA LA BL PI OR RE
JTTA-53NOTC3-2	0.5"	12.7	3"	2	2,000	WH YE GA LA BL PI OR RE
JTTA-106C1-2	0.75"	19.1	1"	1	2,000	WH YE BL
JTTA-106C3-2	0.75"	19.1	3"	1	2,000	WH YE BL
JTTA-115C1-1	1"	25.4	1"	1	1,000	WH YE GA BL PI OR RE
JTTA-115C3-1	1"	25.4	3"	1	1,000	WH YE GA BL PI OR RE



Rectangles + Circle

These thermal-transfer printable cryogenic labels are provided with rectangles and circles to identify the sides and tops of cryo vials and microtubes. They are available in a range of colors and sizes and can be printed with alphanumeric text, as well as 1D and 2D barcodes.

SKU#	Inches (W x H + Dia.)	mm (W x H + Dia.)	Core	Labels Across	Label Sets /Roll	Colors
JTTA-52NOTC1-2	0.94" x 0.5" + 0.437"	23.9 x 12.7 + 11.1	1"	1 set	2,000	WH YE GA BL OR PI LA RE
JTTA-52NOTC3-4	0.94" x 0.5" + 0.437"	23.9 x 12.7 + 11.1	3"	1 set	4,000	WH YE GA BL OR PI LA RE
JTTA-183C1-2WH	1" x 0.5" + 0.375"	25.4 x 12.7 + 9.5	1"	1 set	2,000	WH
JTTA-183C3-4WH	1" x 0.5" + 0.375"	25.4 x 12.7 + 9.5	3"	1 set	4,000	WH
JTTA-70C1-2	1.125" x 0.75" + 0.437"	28.6 x 19.1 + 11.1	1"	1 set	2,000	WH PI
JTTA-70C3-4	1.125" x 0.75" + 0.437"	28.6 x 19.1 + 11.1	3"	1 set	4,000	WH PI
JTTA-161C1-2	1.1875" x 1" + 0.375"	30.2 x 25.4 + 9.5	1"	1 set	2,000	WH YE GA
JTTA-161C3-2	1.1875" x 1" + 0.375"	30.2 x 25.4 + 9.5	3"	1 set	2,000	WH YE GA
JTTA-158C1-2WH	1.25" x 0.5" + 0.437"	31.8 x 12.7 + 11.1	1"	1 set	2,000	WH
JTTA-158C3-4WH	1.25" x 0.5" + 0.437"	31.8 x 12.7 + 11.1	3"	1 set	4,000	WH
JTTA-274C1-2WH	1.25" x 0.5" + 0.375"	31.8 x 12.7 + 9.5	1"	1 set	2,000	WH
JTTA-274C3-4WH	1.25" x 0.5" + 0.375"	31.8 x 12.7 + 9.5	3"	1 set	4,000	WH
JTTA-172C1-2WH	1.25" x 0.625" + 0.4375"	31.8 x 15.9 + 11.1	1"	1 set	2,000	WH
JTTA-172C3-4WH	1.25" x 0.625" + 0.4375"	31.8 x 15.9 + 11.1	3"	1 set	4,000	WH
JTTA-103C1-2WH	1.25" x 0.625" + 0.375"	31.8 x 16 + 9.5	1"	1 set	2,000	WH
JTTA-103C3-4WH	1.25" x 0.625" + 0.375"	31.8 x 16 + 9.5	3"	1 set	4,000	WH



NitroTAG® with SimPEEL™

Patent Pending

NitroTAG labels with SimPEEL technology are even easier to peel and apply. They are provided with an adhesive-free tab that allows the labels to be lifted off their liner with ease, while the glove-friendly design ensures the label won't stick to your hands or gloves during application.



1 After printing SimPEEL™ labels, peel the label from its support liner by lifting the label using the convenient tab.



2 Position the label over the vial or surface while holding the tab.



3 Affix the label then remove the tab as you continue to affix the label.



4 Done.

SKU#	Inches (W x H)	mm (W x H)	Core	Labels Across	Labels/Roll
SPJTТА-29NPC1-1WH	1" x 1"	25.4 x 25.4	1"	1	1,000
SPJTТА-29NPC3-3WH	1" x 1"	25.4 x 25.4	3"	1	3,000
SPJTТА-554NPSBC1-2WH	1.81" x 0.625"	46 x 15.9	1"	1	2,000
SPJTТА-554NPSBC3-2WH	1.81" x 0.625"	46 x 15.9	3"	1	2,000



+1-450-973-9420
cx@ga-international.com
Labtag.com

USA

P.O. Box 2633,
Champlain, NY
12919, United States
1-800-518-0364

CANADA

4455 Rue Louis-B.-Mayer
Laval QC,
H7P 6B5, Canada
1-800-518-0364

NETHERLANDS

Runweg 28A,
5258 BN Berlicum
Netherlands
+31 (0) 73-7370185

UK

Preston Park House, South Road
Brighton, East Sussex
BN1 6SB, United Kingdom
+44 (20) 3769 5683