



-196°C

NitroTAG[®]

Cryogenic Labels for Liquid Nitrogen Storage

Temperature range -196°C to +110°C (-321°F to +230°F)



- Resist temperatures as low as -196°C (-321°F) in liquid or vapor-phase nitrogen
- Strong, permanent cryogenic adhesive allows labels to adhere to wet surfaces
- Withstand spraying and swabbing with alcohol and surface-sanitizing wipes

Thermal-Transfer NitroTAG®

Rectangle

Rectangular NitroTAG® labels are ideal for identifying small-diameter cryovials, microtubes, cardboard or polypropylene cryo/freezer boxes, cell culture plates, and plastic bags. These labels feature a strong, permanent adhesive that remains securely attached, even in ultra-low-temperature environments.

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll	Labels Across	Colors
JTTA-7C1-1	1" x 0.5"	25.4 x 12.7	1"	1,000	1	WH YE GA BL OR RE PI LA SI LE
JTTA-4C1-1WH	1" x 0.75"	25.4 x 19.1	1"	1,000	1	WH
JTTA-174C1-6WH	1" x 1"	25.4 x 25.4	1"	6,000	3	WH
JTTA-29C1-1	1" x 1"	25.4 x 25.4	1"	1,000	1	WH YE GA BL OR RE PI BRG BRB GBL RYE BLA
JTTA-9C1-1	1.25" x 0.875"	31.8 x 22.2	1"	1,000	1	WH YE GA BL OR RE MT
JTTA-58C1-1WH	1.25" x 1"	31.8 x 25.4	1"	1,000	1	WH
JTTA-10C1-1	1.5" x 0.75"	38.1 x 19.1	1"	1,000	1	WH GA OR
JTTA-504C1-1WH	1.5" x 1"	38.1 x 25.4	1"	1,000	1	WH
JTTA-28C1-2	2" x 1"	50.8 x 25.4	1"	2,000	1	WH YE GA BL OR RE PI LA SI BR
JTTA-36C1-1	3" x 1"	76.2 x 25.4	1"	1,000	1	WH YE GA BL OR RE PI LA SI BR
JTTA-100C1-0.5WH	4" x 3"	101.6 x 76.2	1"	500	1	WH

Also available on 3" core for industrial printers.




SimPEEL™

Patented

NitroTAG® labels with SimPEEL™ technology are easy to peel and apply. An adhesive-free tab makes it easy to lift each label from its liner, while the glove-friendly design helps prevent the label from sticking to your hands or gloves.

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll	Colors
SPJTTA-4NPC1-1WH	1" x 0.75"	25.4 x 19.1	1"	1,000	WH
SPJTTA-29NPC1-1WH	1" x 1"	25.4 x 25.4	1"	1,000	
SPJTTA-28NPC1-2WH	2" x 1"	50.8 x 25.4	1"	2,000	

Also available on 3" core for industrial printers.



Rectangle + Circle

These cryogenic labels are provided in rectangles and circles to identify the sides and tops of cryovials and microtubes. Available in a range of colors and sizes, they can be printed with alphanumeric text, 1D barcodes, and 2D barcodes.

SKU#	Inches (W x H + Dia.)	mm (W x H + Dia.)	Core	Labels/Roll	Colors
JTTA-52NOTC1-2	0.94" x 0.5" + 0.437"	23.9 x 12.7 + 11.1	1"	2,000	WH YE GA BL OR RE PI LA
JTTA-183C1-2	1" x 0.5" + 0.375"	25.4 x 12.7 + 9.5	1"	2,000	WH YE GA BL OR LA
JTTA-70C1-2	1.125" x 0.75" + 0.437"	28.6 x 19.1 + 11.1	1"	2,000	WH YE GA BL PI
JTTA-274C1-2WH	1.25" x 0.5" + 0.375"	31.8 x 12.7 + 9.5	1"	2,000	WH
JTTA-172C1-2	1.25" x 0.625" + 0.4375"	31.8 x 15.9 + 11.1	1"	2,000	WH YE GA VI

Also available on 3" core for industrial printers.



SnapPEEL™

© Protected Design

Patent Pending

NitroTAG® cryogenic labels featuring SnapPEEL™ technology are designed to increase both the speed and efficiency of sample labeling and identification. The combined rectangle-and-circle format allows users to quickly label the sides and caps of cryovials and microtubes before storage in liquid nitrogen at -196°C (-321°F) or in ultra-low-temperature freezers at -80°C (-112°F).

SKU#	Inches (W x H + Dia.)	mm (W x H + Dia.)	Core	Labels/Roll	Colors
JTTA-566SBC1-2	0.8" x 0.6" + 0.35"	20.3 x 15.2 + 8.9	1"	2,000	WH YE GA BL OR PI
JTTA-560SBC1-2	0.94" x 0.5" + 0.437"	23.9 x 12.7 + 11.1	1"	2,000	WH YE GA BL OR PI
JTTA-562SBC1-2	1.125" x 0.75" + 0.437"	28.6 x 19.1 + 11.1	1"	2,000	WH YE GA BL OR RE PI

Also available on 3" core for industrial printers.

1

Peel.

2

Stick.

3

Snap.

4

Apply.

Circle

NitroTAG® circle labels are designed to identify vial caps and tube tops. Available in a variety of colors, these cryogenic labels are ideal for color-coding samples and can be printed with text and 2D barcodes for enhanced sample tracking.

SKU#	Inches (Diameter)	mm (Diameter)	Core	Labels/Roll	Colors
JTTA-59NPNOTC1-2	0.25"	6.4	1"	2,000	WH YE GA BL OR PI LA LE
JTTA-61NPNOTC1-4	0.354"	9.0	1"	4,000	WH YE BL OR PI LA SI MT
JTTA-22NOTC1-2	0.438"	11.1	1"	2,000	WH YE GA BL OR RE PI MT
JTTA-53NOTC1-2	0.5"	12.7	1"	2,000	WH YE GA BL OR RE PI LA
JTTA-106C1-2	0.75"	19.1	1"	2,000	WH YE

Also available on 3" core for industrial printers.




Writable NitroTAG®

SnapPEEL™

© Protected Design

Patent Pending

Writable cryo labels with SnapPEEL™ technology, packaged in a convenient, portable dispenser box for easy access, storage, and application. Designed for long-term cryo-preservation under cryogenic laboratory conditions.



ISO 9001:2015 Certified
ISO 22301:2012 Compliant

cx@ga-international.com
Labtag.com