

ParafiTAG™

Labels for Paraffin Wax Blocks

Temperature range -55°C to +70°C (-67°F to +158°F)



SHOP NOW

- Withstand direct exposure to harsh chemicals and solvents
- Adhere to hard-to-stick surfaces, including paraffin wax blocks
- Strong permanent adhesive ensures labels will not wrinkle or peel

ParafiTAG™ thermal-transfer labels are made of durable synthetic material and designed to adhere to paraffin wax blocks as well as other hard-to-label surfaces, such as Teflon, Polytetrafluoroethylene (PTFE), Perfluoroalkoxy (PFA), Parchment paper, and low surface energy plastics. Waterproof, with a strong, extra-permanent adhesive, these labels can withstand exposure to solvents and chemicals, perfectly suited for identifying paraffin wax-embedded tissue cassettes used in histology applications.

Printing with our solvent-resistant ribbons provides a robust chemical-resistant printout capable of withstanding xylene, toluene, alcohols, as well as other chemicals and solvents. These chemical-resistant labels will not fade or peel when immersed in xylene, alcohols, or acetone for up to 15 minutes.

SKU#	Inches (W x H)	mm (W x H)	Core	Labels Across	Labels/Roll
PRF-213C1-5WH	0.75" x 0.2"	19.1 x 5.1	1"	1	5,000
PRF-213C3-10WH	0.75" x 0.2"	19.1 x 5.1	3"	1	10,000
PRF-509C1-4WH	0.75" x 0.25"	19.1 x 6.4	1"	1	4,000
PRF-509C3-8WH	0.75" x 0.25"	19.1 x 6.4	3"	1	8,000
PRF-125C1-4WH	0.9375" x 0.9375"	23.8 x 23.8	1"	4	4,000
PRF-125C3-4WH	0.9375" x 0.9375"	23.8 x 23.8	3"	4	4,000
PRF-67C1-1WH	0.9375" x 0.9375"	23.8 x 23.8	1"	1	1,000
PRF-67C3-1WH	0.9375" x 0.9375"	23.8 x 23.8	3"	1	1,000
PRF-149C1-4WH	1" x 0.276"	25.4 x 7	1"	1	4,000
PRF-149C3-10WH	1" x 0.276"	25.4 x 7	3"	1	10,000
PRF-209C1-4WH	1" x 0.25"	25.4 x 6.4	1"	1	4,000
PRF-209C3-4WH	1" x 0.25"	25.4 x 6.4	3"	1	4,000
PRF-7C1-1WH	1" x 0.5"	25.4 x 12.7	1"	1	1,000
PRF-7C3-3WH	1" x 0.5"	25.4 x 12.7	3"	1	3,000
PRF-242C1-0.5WH	1.75" x 0.75"	44.5 x 19.1	1"	1	500
PRF-242C3-1WH	1.75" x 0.75"	44.5 x 19.1	3"	1	1,000

