

#1 Provider
of Laboratory
identification solutions



Labeling solutions for all lab applications

Cryogenic • Chemical • Autoclave

Major brands of printers provided including

Zebra • cab • DYMO

- Over 20 years of servicing major government, private institutions and research centers.
- Guaranteed compatibility between our labels and printers we provide.
- Learn more with our printer buying guide at Labtag.com/printing-guide.



This table provides an overview of some of the printers that **GA International** offers.

● Standard ○ Option - Not available

	DYMO LabelWriter 450	Zebra TLP 2824 Plus	Zebra ZD620t	cab MACH1	Zebra ZT230	cab SQUIX 4 M
Printing Technology						
• Direct thermal	●	●	●	●	●	●
• Thermal-transfer	-	●	●	●	●	●
Resolution						
• 203 dots per inch	-	●	○	○	○	-
• 300 dots per inch	●	-	●	●	●	●
• 600 dots per inch	-	-	-	-	-	○
Label Dimensions						
• Min. print height	0.5"	0.38"	0.25"	0.16"	0.38"	0.16"
• Max. print width	2.2"	2.24"	4.27"	4.16"	4.09"	4.16"
Capacity						
• Up to 500 labels/day	●	●	●	●	-	-
• More than 500 labels/day	-	-	-	-	●	●
Connectivity						
• USB	●	●	●	●	●	●
• Ethernet	-	○	●	●	○	●
• Wi-Fi	○	-	○	-	○	●
• Bluetooth	-	-	●	-	-	-
Add-ons						
• Peeler	-	○	○	○	○	○
• Cutter	-	○	○	○	○	○
• Rewinder	-	-	-	-	-	○
• External unwinder	-	○	-	○	-	○
Applications						
• Cryogenic labels (no alcohol resistance)	●	●	●	●	●	●
• Cryogenic labels (alcohol resistance)	●	●	●	●	●	●
• Cryogenic labels (frozen surfaces)	●	●	●	●	●	●
• Autoclave labels	-	●	●	●	●	●
• Chemical-resistant (solvents, stains)	-	●	●	●	●	●
• Most laboratory applications	●	●	●	●	●	●



Let our scientific experts help you choose the best printer and labels for your labeling needs and provide you with full after-sales support.