



Featured: GA Cody-3

Thermal Printers for the Lab



Labeling solutions for all lab applications

Cryogenic • Chemical • Autoclave

Major brands of printers provided including

Zebra • cab • GA International

- › Over 20 years of servicing major government, private institutions and research centers.
- › Guaranteed compatibility between our labels and printers we provide.
- › Learn more with our printer buying guide at Labtag.com/printing-guide.



Featured: GA DiTi 203iD

This table provides an overview of some of the printers that **GA International** offers.

● Standard ○ Option - Not available

	GA Cody-3	GA DiTi 203iD	Zebra ZD621t	cab MACH1	Zebra TLP 2824 Plus	cab SQUIX 4 M	Zebra ZT230
Printing Technology							
• Direct thermal	●	●	●	●	●	●	●
• Thermal-transfer	●	-	●	●	●	●	●
Resolution							
• 203 dots per inch	-	●	○	○	●	-	○
• 300 dots per inch	●	-	●	●	-	●	●
• 600 dots per inch	-	-	-	-	-	○	-
Label Dimensions							
• Min. print height	0.25"	0.5"	0.25"	0.16"	0.38"	0.16"	0.38"
• Max. print width	4.16"	3.35"	4.27"	4.16"	2.24"	4.16"	4.09"
Capacity							
• More than 500 labels/day	-	-	-	-	-	●	●
Connectivity							
• USB	●	●	●	●	●	●	●
• Ethernet	●	-	●	●	○	●	○
• Wi-Fi	-	-	-	-	-	●	○
• Bluetooth	-	-	●	-	-	●	-
Add-ons							
• Peeler	-	-	○	-	○	○	○
• Cutter	-	-	-	○	○	○	○
• Rewinder	-	-	-	-	-	○	-
• External unwinder	-	-	○	○	○	○	○
Applications							
• Cryogenic labels	●	●	●	●	●	●	●
• Cryogenic labels (frozen surfaces)	●	●	●	●	●	●	●
• Autoclave labels	●	-	●	●	●	●	●
• Chemical-resistant (solvents, stains)	●	●	●	●	●	●	●
• Most laboratory applications	●	●	●	●	●	●	●



*Please Note, printers are not shown to scale

Let our scientific experts help you choose the best printer and labels for your labeling needs and provide you with full after-sales support.