

Steri-ThermoTAG™

Autoclave-Resistant Thermal-Transfer Labels

Temperature range from -60°C to +150°C (-76°F to +302°F)



- Able to withstand a wide temperature range (-60°C/-76°F to +150°C/+302°F)
- Designed for use in steam autoclaves, dry heat ovens, ethylene oxide sterilizers, and gamma irradiators
- Versatile, can adhere to a variety of surfaces, including metal, plastic, and glass
- Ideal for labeling surgical instruments, plasticware, glassware, and blister packs
- Compatible with all common thermal-transfer printer brands

Recommended Ribbons: Smudge & Alcohol Proof RR-Class

Shop on www.Labtag.com/steri-thermotag



Steri-ThermoTAG™ Autoclave-Resistant Thermal-Transfer Labels

Permanent

These permanent autoclavable thermal-transfer labels have been tested to endure up to 20 autoclave cycles, perfectly suited for the identification of containers and instruments that will undergo multiple sterilization cycles.

SKU# (White)	Inches (W x H)	mm (W x H)	Core	Labels/Roll
AUTT-28C1-2WH	2" x 1"	50.8 x 25.4	1"	2,000
AUTT-28C3-2WH	2" x 1"	50.8 x 25.4	3"	2,000
AUTT-200C1-1WH	2.5" x 1.25"	63.5 x 31.8	1"	1,000
AUTT-200C3-2WH	2.5" x 1.25"	63.5 x 31.8	3"	2,000



Removable

Coated with a glove-friendly removable adhesive, these thermal-transfer steam autoclave labels can be readily removed without leaving any glue residue behind, ideal for labeling labware, stainless steel surfaces, and containers for sterilization.

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll
AUTR-9C1-1WH	1.25" x 0.875"	31.8 x 22.2	1"	1,000
AUTR-9C3-3WH	1.25" x 0.875"	31.8 x 22.2	3"	3,000
AUTR-28C1-1WH	2" x 1"	50.8 x 25.4	1"	1,000
AUTR-28C3-2WH	2" x 1"	50.8 x 25.4	3"	2,000
AUTR-56C1-1WH	3" x 2"	76.2 x 50.8	1"	1,000
AUTR-56C3-1WH	3" x 2"	76.2 x 50.8	3"	1,000



Blackout

Blackout Steri-ThermoTAG™ labels are ideal for concealing sensitive information on containers that need to be autoclaved before disposal.

SKU#	Inches (W x H)	mm (W x H)	Core	Labels/Roll
BOA-19C1-1SI	1.3" x 0.625"	33 x 15.9	1"	1,000
BOA-19C3-3SI	1.3" x 0.625"	33 x 15.9	3"	3,000
BOA-28C1-1SI	2" x 1"	50.8 x 25.4	1"	1,000
BOA-28C3-2SI	2" x 1"	50.8 x 25.4	3"	2,000



Steri-ThermoTAG is also available in color options for some sizes, please visit Labtag.com/steri-thermotag
For custom color and sizing options, please contact: support@ga-international.com