

UbiCling™

Static Cling Labels & Tapes



Labels

Unique specialty labels that swiftly adhere to metal, plastic, cardboard, wooden, and glass surfaces using static-cling like properties. UbiCling labels are easily repositionable and will not leave behind any residue when removed. Thermal-transfer printable and writable with markers and pens, they are ideal for identifying most non-frosted glass surfaces such as labware, Pyrex bottles, windows, aquariums, as well as other smooth surfaces, including bags, coated paper and carton, stainless steel, laminated wood products, Teflon, epoxy resin, and polished carbon fiber.

White & Colors

SKU#	Inches	mm	Core	Labels/Roll	Colors Available
GFS-118NOTC1-1WH	1" x 1.25"	25.4 x 31.8	1"	1,000	(WH)
GFS-118NOTC3-1WH	1" x 1.25"	25.4 x 31.8	3"	1,000	(WH)
GFS-551NOTC1-0.5	2.125" x 1.375"	54 x 35	1"	500	(WH) (YE) (BL) (PI) (OR) (GA) (RE)
GFS-551NOTC3-0.5	2.125" x 1.375"	54 x 35	3"	500	(WH) (YE) (BL) (PI) (OR) (GA) (RE)



Clear

SKU#	Inches	mm	Core	Labels/Roll
GFC-53NOTC1-2	0.5" Circle	12.7	1"	2,000
GFC-53NOTC3-4	0.5" Circle	12.7	3"	4,000
GFC-27NOTC1-0.5	2" x 2"	50.8 x 50.8	1"	500
GFC-27NOTC3-1	2" x 2"	50.8 x 50.8	3"	1,000
GFC-28NOTC1-1	2" x 1"	50.8 x 25.4	1"	1,000
GFC-28NOTC3-2	2" x 1"	50.8 x 25.4	3"	2,000






Tapes

Our static-cling tape is provided with GatorCUT™ technology that makes it fully portable, able to be hand-torn to the desired length effortlessly. Its glove-friendly, repositionable static-cling design makes it perfect for labeling various surfaces, most non-frosted glass surfaces, galvanized, epoxy coated, painted, and stainless steel, as well as polished, ceramic, porcelain surfaces.


White

SKU#	Inches x Feet	Metric	Core
TGS-25CP1-50WH	1" x 50'	25.4mm x 15m	1"
TGS-25CP3-50WH	1" x 50'	25.4mm x 15m	3"



Clear

SKU#	Inches x Feet	Metric	Core
TGC-25CP1-50	1" x 50'	25.4mm x 15m	1"
TGC-25CP3-50	1" x 50'	25.4mm x 15m	3"



 [Shop on Labtag.com/ubicling](https://www.Labtag.com/ubicling)