


Vial Cap Sizing Guide

Use our handy Vial Cap Sizing Guide to ensure the best fit and allow optimal performance of the labels. Sizes of caps may vary according to the manufacturer. Users are encouraged to assess the suitability of sizes and materials according to their specific needs. Please visit www.Labtag.com to obtain expert advice and how to obtain free samples.

A	0.2 ml flat-top PCR tubes
B	0.5 ml flat-top microtubes
C	1.5 ml flat-top microtubes
D	1-2 ml flat-top cryovials
E	1-5 ml indented-cap cryovial
F	15 ml flat top sample tubes
G	50 ml flat top sample tubes



Thermal-Transfer Labels

Label	Label Size	Recommended Tube Vial Type	PikaTAG® Compatible
JTTA-294	0.236" x 0.236" / 6mm x 6mm	B	No*
JTTA-59	0.25" / 6.4mm	A,B	No*
JTTA-288	0.25" / 6.4mm	A,B	No*
JTTA-61	0.35" / 9mm	B	No*
JTTA-255	0.375" x 0.375" / 9.5mm x 9.5mm	B	No*
JTTA-22	0.438" / 11.1mm	C	No
JTTA-204	0.47" x 0.47" / 11.9mm x 11.9mm	C	No
JTTA-53	0.5" / 13mm	F	No
JTTA-293	0.5" / 13mm	F	No
JTTA-500	0.5" / 13mm	F	No
JTTA-230	0.6" x 0.6" + 0.375" / 15.2mm x 15.2mm + 9.5mm	C,D	No*
JTTA-208	0.7885" x 1.375" + 0.4375" / 20mm x 35 mm + 11.1mm	C,D	No
JTTA-52	0.94" x 0.5" + 0.437" / 23.9mm x 12.7mm + 11.1mm	C,D	No
JTTA-228	1" x 0.375" + 0.375" / 25.4mm x 9.5mm + 9.5mm	C,D	No*
JTTA-183	1" x 0.5" + 0.375" / 25.4mm x 12.7mm + 9.5mm	C,D	No*
JTTA-170	1.125" x 0.75" + 0.437" / 28.6mm x 19.1mm + 11.1mm	C,D	No
JTTA-70	1.125" x 0.75" + 0.437" / 28.6mm x 19.1mm + 11.1mm	C,D	No
JTTA-124	1.18" x 0.28" + 0.37" / 30mm x 7.1mm + 9.4mm	C,D	No*
JTTA-161	1.1875" x 1" + 0.375" / 30.2mm x 25.4mm + 9.5mm	C,D	No*
JTTA-171	1.25" x 0.625" + 0.4375" / 31.8mm x 16mm + 11.1mm	C,D	No
JTTA-266	1.25" x 0.625" + 0.375" / 31.8mm x 15.9mm + 9.5mm	C,D	No*
JTTA-172	1.25" x 0.625" + 0.4375" / 31.8mm x 16mm + 11.1mm	C,D	No
JTTA-274	1.25" x 0.5" + 0.375" / 31.8mm x 12.7mm + 9.5mm	C,D	No*
JTTA-103	1.25" x 0.625" + 0.375" / 31.8mm x 16mm + 9.5mm	C,D	No*
JTTA-158	1.25" x 0.5" + 0.437" / 31.8mm x 12.7mm + 11.1mm	C,D	No
JTTA-197	1.438" x 0.55" + 0.375" / 36.5mm x 14mm + 9.5mm	C,D	No*
JTTA-514	1.73" x 0.63" + 0.374" / 44mm x 16mm + 9.5mm	C,D	No
JTTA-292	0.98" x 0.59" + 0.33" / 24.9mm x 15mm + 8.4mm	C	No*
GANA-52	0.94" x 0.5" + 0.437" / 23.9mm x 12.7mm + 11.1mm	C,D	No
HBCL-61	0.354" / 9mm	C,D	No
VCL-61NPNOT	0.354" / 9mm	E	Yes
VCSL-555NOT	1.125" x 0.625" + 0.354" / 28.6mm x 15.9mm + 9mm	E	Yes
FJT-290	2" x 1" + 0.75" / 50.8mm x 25.4mm + 19.1mm	F, G	No
FJT-556	2" x 1" + 1" / 50.8mm x 25.4mm + 25.4mm	G	No
UC-52	0.94" x 0.5" + 0.437" / 23.9mm x 12.7mm + 11.1mm	C, D	No

*A deep-freeze (-80°C) version of these labels (FJT class) is available which is compatible with PikaTAG®



Direct Thermal Labels

Label	Label Size	Recommended Tube Vial Type	PikaTAG® Compatible
DFP-161	1.1875" x 1" + 0.375" / 30.2mm x 25.4mm + 9.5mm	C,D	No
DFP-103	1.25" x 0.625" + 0.375" / 31.8mm x 16mm + 9.5mm	C,D	No
DFP-204	0.47" x 0.47" / 11.9mm x 11.9mm	C,D	No
DFP-52	0.94" x 0.5" + 0.437" / 23.9mm x 12.7mm + 11.1mm	C,D	No
DFP-61	0.354" / 9mm	C,D	Yes
DTU-61	0.35" / 9mm	C, D	No
DTU-52	0.94" x 0.5" + 0.437" / 23.9mm x 12.7mm + 11.1mm	C, D	No



DYMO-Compatible Labels

Label	Label Size	Recommended Tube Vial Type	PikaTAG® Compatible
ED1F/EF1F-030	0.35" / 9mm	C,D	Yes
ED1F/EF1F-050	0.5" / 13mm	F	No
EC1F/EG1F-072	0.5" x 1.04" + 0.375" / 12.7mm x 26mm + 9.5mm	C,D	Yes
EC1F/EG1F-081	1.04" x 0.5" + 0.375" / 26.4mm x 12.7mm + 9.5mm	C,D	Yes
ED1F-079	0.79" x 0.2" + 0.25" / 20.1mm x 5.1mm + 6.4mm	A	Yes
EDCS-030	0.35" / 9mm	B	No
EDCS-072	0.5" / 13mm	F	No

Inkjet Labels

Label	Label Size	Recommended Tube Vial Type	PikaTAG® Compatible
CIJ-52	0.94" x 0.5" + 0.437" / 23.9mm x 12.7mm + 11.1mm	C, D	Yes
CIJ-183	1.0" x 0.5" + 0.375" / 25.4mm x 12.7mm + 9.5mm	C, D	Yes
CIJS-10	0.35" / 9mm	C, D	Yes
CIJS-48	0.43" / 11mm	F	Yes

Laser Printer Labels

Label	Label Size	Recommended Tube Vial Type	PikaTAG® Compatible
CL-48	0.43" / 11mm	C,D	No
A4CL-48	0.43" / 11mm	C,D	No



SnapPEEL Labels



Label	Label Size	Recommended Tube Vial Type	PikaTAG® Compatible
JTTA-560	0.94" x 0.5" + 0.437" / 23.9mm x 12.7mm + 11.1mm	C, D	No
JTTA-562	1.125" x 0.75" + 0.437" / 28.6mm x 19.1mm + 11.1mm	C, D	No
FJT-564	0.75" x 0.25" + 0.25" / 19.1mm x 6.4mm + 6.4mm	A, B	No
DTU-560	0.94" x 0.5" + 0.437" / 23.9mm x 12.7mm + 11.1mm	C, D	No
LCS-125	0.94" x 0.5" + 0.437" / 23.9mm x 12.7mm + 11.1mm	C, D	No

Writable Only Labels

Label	Label Size	Recommended Tube Vial Type	PikaTAG® Compatible
LT-9	0.35" / 9mm	B	Yes
LT-11	0.44" / 11mm	C	No
LT-13	0.5" / 13mm	F	No
LT-19	0.75" / 19mm	G	No
WRPLT-7	0.28" / 7mm	A	Yes
JTA-13X6A	0.5" / 13mm	F	No
JTRA-9	0.35" / 9mm	B	Yes
JTRA-11	0.44" / 11mm	C	No
JTRA-13	0.5" / 13mm	F	No



+1-450-973-9420
 cx@ga-international.com
 Labtag.com

USA

P.O. Box 2633,
 Champlain, NY
 12919, United States
 1-800-518-0364

CANADA

4455 Rue Louis-B.-Mayer,
 Laval QC,
 H7P 6B5, Canada
 1-800-518-0364

NETHERLANDS

Runweg 28A,
 5258 BN Berlicum
 Netherlands
 +31 (0) 73-7370185

UK

Preston Park House, South Road
 Brighton, East Sussex
 BN1 6SB, United Kingdom
 +44 (20) 3769 5683