

New!

# Zesti-DTermo™

## DYMO-Compatible Heat-Resistant Labels

- Compatible with DYMO LabelWriter® printers such as the 450, 450 Turbo and 4XL
- Ideal for identifying vials, tubes, microtubes, racks, and other laboratory containers
- Available in various sizes, offering resistance to high temperatures as well as incidental exposure to some common chemicals

Patent Pending



Our DYMO-compatible heat-resistant labels withstand elevated temperatures, up to +95°C (+203°F), and are well suited for heat inactivation and high heat treatment protocols. In contrast to standard DYMO-compatible labels, these labels will not blacken from exposure to high heat and will remain legible and scannable with barcode scanners. Their strong permanent adhesive allows them to adhere to metal, glass, and plastic surfaces.



Temperature range from -80°C to +95°C (-112°F to +203°F)

SKU#	Inches (W x H)	mm (W x H)	Labels Across	Labels/Roll
EDAP-040WH	0.5" x 1"	12.7 x 25.4	2	1,000
EDAP-070WH	1" x 1"	25.4 x 25.4	1	750
EDAP-082WH	2.25" x 1.25"	57 x 32	1	500

