


Vial Cap Sizing Guide

Use our handy Vial Cap Sizing Guide to ensure the best fit and allow optimal performance of the labels. Sizes of caps may vary according to the manufacturer. Users are encouraged to assess the suitability of sizes and materials according to their specific needs.

Please visit www.Labtag.com to obtain expert advice and how to obtain free samples.

A	0.2 ml flat-top PCR tubes
B	0.5 ml flat-top microtubes
C	1.5 ml flat-top microtubes
D	1-2 ml flat-top cryovials
E	1-5 ml indented-cap cryovial
F	15 ml flat top sample tubes
G	50 ml flat top sample tubes



Writable Only Labels

Label	Label Size	Recommended Tube Vial Type	PikaTAG® Compatible
LT-9	0.35" / 9mm	B	Yes
LT-11	0.44" / 11mm	C	No
LT-13	0.5" / 13mm	F	No
LT-19	0.75" / 19mm	G	No
WRPLT-7	0.28" / 7mm	A	Yes
JTA-13X6A	0.5" / 13mm	F	No
JTRA-9	0.35" / 9mm	B	Yes
JTRA-11	0.44" / 11mm	C	No
JTRA-13	0.5" / 13mm	F	No

Thermal-Transfer Labels

Label	Label Size	Recommended Tube Vial Type	PikaTAG® Compatible
JTTA-294	0.236" x 0.236" / 6mm x 6mm	B	No*
JTTA-59	0.25" / 6.4mm	A,B	No*
JTTA-288	0.25" / 6.4mm	A,B	No*
JTTA-61	0.35" / 9mm	B	No*
JTTA-255	0.375" x 0.375" / 9.5mm x 9.5mm	B	No*
JTTA-22	0.438" / 11.1mm	C	No
JTTA-204	0.47" x 0.47" / 11.9mm x 11.9mm	C	No
JTTA-53	0.5" / 13mm	F	No
JTTA-293	0.5" / 13mm	F	No
JTTA-500	0.5" / 13mm	F	No
JTTA-230	0.6" x 0.6" + 0.375" / 15.2mm x 15.2mm + 9.5mm	C,D	No*
JTTA-208	0.7885" x 1.375" + 0.4375" / 20mm x 35 mm + 11.1mm	C,D	No
JTTA-52	0.94" x 0.5" + 0.437" / 23.9mm x 12.7mm + 11.1mm	C,D	No
JTTA-228	1" x 0.375" + 0.375" / 25.4mm x 9.5mm + 9.5mm	C,D	No*
JTTA-183	1" x 0.5" + 0.375" / 25.4mm x 12.7mm + 9.5mm	C,D	No*
JTTA-170	1.125" x 0.75" + 0.437" / 28.6mm x 19.1mm + 11.1mm	C,D	No
JTTA-70	1.125" x 0.75" + 0.437" / 28.6mm x 19.1mm + 11.1mm	C,D	No
JTTA-124	1.18" x 0.28" + 0.37" / 30mm x 7.1mm + 9.4mm	C,D	No*
JTTA-161	1.1875" x 1" + 0.375" / 30.2mm x 25.4mm + 9.5mm	C,D	No*
JTTA-171	1.25" x 0.625" + 0.4375" / 31.8mm x 16mm + 11.1mm	C,D	No
JTTA-266	1.25" x 0.625" + 0.375" / 31.8mm x 15.9mm + 9.5mm	C,D	No*
JTTA-172	1.25" x 0.625" + 0.4375" / 31.8mm x 16mm + 11.1mm	C,D	No
JTTA-274	1.25" x 0.5" + 0.375" / 31.8mm x 12.7mm + 9.5mm	C,D	No*
JTTA-103	1.25" x 0.625" + 0.375" / 31.8mm x 16mm + 9.5mm	C,D	No*
JTTA-158	1.25" x 0.5" + 0.437" / 31.8mm x 12.7mm + 11.1mm	C,D	No
JTTA-197	1.438" x 0.55" + 0.375" / 36.5mm x 14mm + 9.5mm	C,D	No*
JTTA-514	1.73" x 0.63" + 0.374" / 44mm x 16mm + 9.5mm	C,D	No
JTTA-292	0.98" x 0.59" + 0.33" / 24.9mm x 15mm + 8.4mm	C	No*
GANA-52	0.94" x 0.5" + 0.437" / 23.9mm x 12.7mm + 11.1mm	C,D	No
HBCL-61	0.354" / 9mm	C,D	No
VCL-61NPNOT	0.354" / 9mm	E	Yes
VCSL-555NOT	1.125" x 0.625" + 0.354" / 28.6mm x 15.9mm + 9mm	E	Yes

* Deep-freeze (-80°C) version of these labels are available that are compatible with PikaTAG

DYMO-Compatible Labels

Label	Label Size	Recommended Tube Vial Type	PikaTAG® Compatible
ED1F/EF1F-030	0.35" / 9mm	C,D	Yes
ED1F/EF1F-050	0.5" / 13mm	F	No
EC1F/EG1F-072	0.5" x 1.04" + 0.375" / 12.7mm x 26mm + 9.5mm	C,D	Yes
EC1F/EG1F-081	1.04" x 0.5" + 0.375" / 26.4mm x 12.7mm + 9.5mm	C,D	Yes
ED1F-079	0.79" x 0.2" + 0.25" / 20.1mm x 5.1mm + 6.4mm	A	Yes
EDCC-030	0.35" / 9mm	B	Yes
EDCC-072	0.5" / 13mm	F	Yes

Inkjet Labels

Label	Label Size	Recommended Tube Vial Type	PikaTAG® Compatible
ADA-52	0.94" x 0.5" + 0.437" / 23.9mm x 12.7mm + 11.1mm	C,D	No

Laser Printer Labels

Label	Label Size	Recommended Tube Vial Type	PikaTAG® Compatible
CL-48	0.43" / 11mm	C,D	No
A4CL-48	0.43" / 11mm	C,D	No

Direct Thermal Labels

Label	Label Size	Recommended Tube Vial Type	PikaTAG® Compatible
DFP-161	1.1875" x 1" + 0.375" / 30.2mm x 25.4mm + 9.5mm	C,D	No
DFP-103	1.25" x 0.625" + 0.375" / 31.8mm x 16mm + 9.5mm	C,D	No
DFP-204	0.47" x 0.47" / 11.9mm x 11.9mm	C,D	No
DFP-52	0.94" x 0.5" + 0.437" / 23.9mm x 12.7mm + 11.1mm	C,D	No
DFP-61	0.354" / 9mm	C,D	Yes

Labtag.com/vial-cap-labels